



PRE-BID DOCUMENT

| | |
|-------------------------------|--|
| PRE-BID NUMBER: | ZNB 6745/2024-H |
| DESCRIPTION: | THE SUPPLY AND DELIVERY OF HOSPITAL SUNDRIES FOR VARIOUS INSTITUTIONS |
| PERIOD | THREE (03) YEARS CONTRACT |
| PRE-BID MEETING DATE: | 20 March 2024 |
| PRE-BID MEETING TIME: | 10H00 -13H00 |
| PRE-BID MEETING VENUE: | <p>Microsoft Teams meeting</p> <p>Join on your computer, mobile app or room device</p> <p>Click here to join the meeting</p> <p>Meeting ID: 357 229 517 920</p> <p>Passcode: LP3JQs</p> <p>Download Teams Join on the web</p> <p>Or call in (audio only)</p> <p>+27 21 834 5669,,791557038# South Africa, Cape Town</p> <p>Phone Conference ID: 791 557 038#</p> |

Bidders Notes :

Prospective bidders must read this document and make prior notes for further discussions at the scheduled pre-bid specification meeting.

SPECIFICATIONS THE SUPPLY AND DELIVERY OF HOSPITAL SUNDRIES FOR VARIOUS INSTITUTIONS

SPECIFICATIONS LIST FOR HOSPITAL SUNDRIES

| GENERAL SUNDRIES | | | | |
|--------------------------------------|------------|-------------|--------------|--|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| PERSONAL IDENTIFICATION BANDS | 1 | | 30 028 07 | Personal Identification Band - Adult Yellow Locking device from 12cm - 25cm |
| | 2 | | 30 028 08 | Personal Identification Band - Adult Pink Locking device from 12cm - 25cm |
| | 3 | | 30 028 09 | Personal Identification Band - Adult White Locking device from 12cm - 25cm |
| | 4 | | 30 028 10 | Personal Identification Band - Child White Locking device from 9cm - 21cm |
| | 5 | | 30 028 12 | Personal Identification Band - Child Pink Locking device from 9cm - 21cm |
| | 6 | | 30 028 13 | Personal Identification Band - Child Yellow Locking device from 9cm - 21cm |
| | 7 | | 30 028 11 | Personal Identification Band - Neonate and Premature Babies White Locking device from 6cm - 12cm |
| DISPOSABLE UNDERWEAR | 8 | | 39 054 13 | Disposable underwear - Small Briefs Waist circumference: 80-90cm Leg circumference: 28-60cm Waist-crotch: ≥ 30cm |
| | 9 | | 39 054 14 | Disposable underwear - Medium Briefs Waist circumference: 90-100cm Leg circumference: 32-64cm Waist-crotch: ≥ 32cm |
| | 10 | | 39 054 15 | Disposable underwear - Large Briefs Waist circumference: 100-120cm Leg circumference: 36-72cm Waist-crotch: ≥ 35cm |
| | 11 | | 39 054 16 | Disposable underwear - X-Large Briefs Waist circumference: 120-140cm Leg circumference: 36-76cm Waist-crotch: ≥ 40cm |
| | 12 | | 39 054 17 | Disposable underwear - XX-Large Briefs Waist circumference: 130cm-150cm Leg circumference: 40cm-80cm Waist-crotch: 44cm |
| | 13 | | 39 054 18 | Disposable underwear - XXX-Large Briefs Waist circumference: 140cm-160cm Leg circumference: 44cm-88cm Waist-crotch: 46cm |
| FLUSHABLE COVERS | 14 | | 39 050 11 | Bedpan Cover - Flushable Size: 30.5cm x 7.6cm x 40.6cm |
| | 15 | | 39 050 12 | Urinal Covers - Flushable Size: 15.22cm x 7.6cm x 30.5cm |
| PROTECTIVE COVERS | 16 | | USPcover | Ultrasound probe cover Size: ± 10cm x ± 60cm |
| | 17 | | 30 461 10 | Temperature probe cover Length: 15-20cm |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|-----------------|-----|------|-----------|---|
| WIPES | 18 | | 30 568 79 | Alcohol Swabs - 24mm x 30mm |
| | 19 | | AdhWipe | Wipe to remove adhesiveness from the skin. Size: 40cm x 40cm |
| GELS & CREAMS | 20 | | 30 344 30 | Ultrasound Gel with dispenser 200ml |
| | 21 | | 30 344 32 | Ultrasound Gel Jerry Can 5 Litre |
| | 22 | | LubGel | Lubricating (KY) Jelly 2.5g Sachets |
| | 23 | | 53 026 10 | Anaesthetic Gel 2% 20ml Tube |
| ECG ELECTRODES | 24 | | 30 248 03 | Solid Gel ECG Electrode - Adult |
| | 25 | | 30 248 04 | Solid Gel ECG Electrode - Paediatric < 50mm wide |
| | 26 | | 30 248 09 | Sponge ECG Electrode - Paediatric < 50mm wide |
| | 27 | | 30 248 18 | Sponge ECG Electrode - Adult |
| | 28 | | 30 248 27 | Tape ECG Electrode - Paediatric < 50mm wide |
| | 29 | | 30 248 37 | Tape ECG Electrode - Neonate < 30mm wide |
| | 30 | | 30 248 38 | Tape ECG Electrode - Adult |
| PLASTIC APRONS | 31 | | 39 010 15 | Full Body Apron Medium - General use Length from neck: ± 110cm Colour: Blue |
| | 32 | | 39 010 16 | Full Body Apron Large - General use Length from neck: ± 130cm Colour: Blue |
| | 33 | | 39 010 18 | Full Body Apron Medium - Sluice room Length from neck: ± 110cm Colour: Red |
| | 34 | | 39 010 19 | Full Body Apron Large - Sluice room Length from neck: ± 130cm Colour: Red |
| | 35 | | 39 010 21 | Full Body Apron Medium - Infectious areas Length from neck: ± 110cm Colour: Yellow |
| | 36 | | 39 010 22 | Full Body Apron Large - Infectious areas Length from neck: ± 130cm Colour: Yellow |
| | 37 | | 39 010 24 | Full Body Apron Medium - Clean areas (kitchen) Length from neck: ± 110cm Colour: White |
| | 38 | | 39 010 25 | Full Body Apron Large - Clean areas (kitchen) Length from neck: ± 130cm Colour: White |
| | 39 | | 39 010 26 | Surgical Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Light Green |
| | 40 | | 54 020 10 | Laundry Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Red |
| NON-DISPOSABLES | 41 | | 30 542 06 | Stethoscope – Infant Diaphragm: ± 2.7cm |
| | 42 | | 30 542 07 | Stethoscope – Paediatric Diaphragm: ± 3.3cm |
| | 43 | | 30 542 08 | Stethoscope – Adult Diaphragm: ± 5cm |
| | 44 | | 30 244 00 | Thermometer – Digital temperature ranges of at least: 34°C to 42°C. |
| | 45 | | ND-PB | Pressure Bag – Non-Disposable |
| | 46 | | ET-M | Manometer - Endotracheal and Supraglottic Cuff Pressure Monitor |
| PASTEURISATION | 47 | | 30 047 00 | Pasteurisation Bottle - Glass. Volume: 125ml |
| | 48 | | 30 047 01 | Lid for Glass Pasteurisation Bottle 48mm single thread |
| | 49 | | 30 047 02 | Pasteurisation and Human Milk Storage Bottle - BPA-free Plastic Volume: 130ml ; Neck Diameter: 50mm |
| | 50 | | 30 047 03 | Lids with foil seal for Human Milk Storage Bottles - Reusable |
| | 51 | | 30 052 09 | Bottle Brush: Length including brush: 40cm. Brush: 10cm, Handle: 30cm Flexible bristles: 5cm |
| | 52 | | PL1 | Small-scale Pasteuriser Labels Size: 28mm x 89mm |
| | 53 | | PL2 | Plasticised Pasteuriser Labels Size: 100mm x 25mm |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|--|-----|------|-----------|--|
| NEONATAL | 54 | | 30 445 23 | Litmus Paper Dark Blue |
| | 55 | | 30 445 24 | Litmus Paper Red |
| | 56 | | 30 445 25 | pH Indicator Test strips |
| | 57 | | 30 460 51 | Phototherapy Eye-Shield - Preemie , head circumference ± 26-32cm |
| | 58 | | 30 460 52 | Phototherapy Eye-Shield - Regular , head circumference ± 33-38 cm |
| | 59 | | 30 072 01 | Umbilical Cord Clamp. ± 5cm long |
| | 60 | | 30 326 10 | Amniotomy hook - 25cm |
| CONNECTORS | 61 | | 30 089 23 | Y-Shaped Tubing Connector - Large OD: ±10mm All sides must be of equal length – 30mm length from point to junction |
| | 62 | | 30 089 36 | Straight Tubing Connector - Open-ended . Double tapered Small OD: 7-11mm; ID: 5mm; Total length: 6cm |
| | 63 | | 30 089 37 | Straight Tubing Connector - Open-ended . Double tapered Large OD: 10 - 14mm ID: 8mm Total length: 6cm |
| | 64 | | 30 534 80 | Spigot - Small Diameter: 9mm ; Length: ±40mm |
| | 65 | | 30 534 82 | Spigot - Large Diameter: 10mm ; Length: ± 52mm |
| DRESSING PACKS | 66 | | 30 432 15 | Large Tray Dressing Pack for Minor Surgical Procedures |
| | 67 | | 30 432 16 | Multi-Pack Tray Dressing Pack for Minor Surgical Procedures |
| | 68 | | 30 432 17 | Catheter Tray Dressing Pack |
| | 69 | | 30 432 18 | Maternity Swabbing Tray Dressing Pack |
| COTTON WOOL AND CLEANSING SWABS | 70 | | 30 645 05 | Cotton wool roll – 1st Grade 50g roll |
| | 71 | | 30 645 06 | Cotton wool roll – 2nd Grade 50g roll |
| | 72 | | 30 568 03 | Cotton Wool Absorbent Balls – Sterile 5 balls Sterile in an easy to open peel pouch |
| | 73 | | 30 568 04 | Cotton Wool Absorbent Balls – Unsterile Box of 1000 balls |
| | 74 | | 30 568 05 | Cleansing ball swabs - X-Ray detectable, Unsterile , Diameter: 65mm |
| | 75 | | 30 568 06 | Cleansing ball swabs - X-Ray detectable, Sterile ; Diameter: 65mm |

| PHARMACY | | | | |
|------------------|-----|------|-----------|--|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| CAPSULE BAGS | 1 | | 74 010 05 | Re-sealable Tablets/Capsules Bags Size: 1 Length: ±100mm Width: ±104mm Lip: ±10mm Overprinted conventional format |
| | 2 | | 74 010 15 | Re-sealable Tablets/Capsules Bags Size: 2 Length: ±120mm Width: ±104mm Lip: ±10mm Overprinted conventional format |
| | 3 | | 74 010 45 | Re-sealable Tablets/Capsules Bags Size: 3 Length: ±160mm Width: ±104mm Lip: ±10mm Overprinted conventional format |
| TABLET VIALS | 4 | | 74 200 11 | Tablet Vial with hinged closure 10ml |
| | 5 | | 74 200 21 | Tablet Vial with hinged closure 20ml |
| | 6 | | 74 200 31 | Tablet Vial with hinged closure 30ml |
| | 7 | | 74 200 51 | Tablet Vial with hinged closure 50ml |
| | 8 | | 74 200 61 | Tablet Vial with hinged closure 75ml |
| | 9 | | 74 200 71 | Tablet Vial with hinged closure 100ml |
| BROWN PAPER BAGS | 10 | | 40 017 40 | Bags Brown paper Size: 1/2 Width: 75mm ; Gusset: 45mm ; Length: 150mm |
| | 11 | | 40 017 42 | Bags Brown paper Size: 1 Width: 80mm ; Gusset: 60mm ; Length: 180mm |
| | 12 | | 40 017 46 | Bags Brown paper Size: 2 Width: 95mm ; Gusset: 60mm ; Length: 110mm |
| | 13 | | 40 017 48 | Bags Brown paper Size: 3 Width: 128mm ; Gusset: 74mm ; Length: 195mm |
| | 14 | | 40 017 50 | Bags Brown paper Size: 4 Width: 124mm ; Gusset: 74mm ; Length: 225mm |
| | 15 | | 40 017 58 | Bags Brown paper Size: 6 Width: 155mm ; Gusset: 90mm ; Length: 301mm |
| | 16 | | 40 017 60 | Bags Brown paper Size: 8 Width: 165mm ; Gusset: 102mm ; Length: 301mm |
| | 17 | | 40 017 62 | Bags Brown paper Size: 10 Width: 165mm ; Gusset: 102mm ; Length: 345mm |
| | 18 | | 40 017 66 | Bags Brown paper Size: 12 Width: 180mm ; Gusset: 102mm ; Length: 335mm |
| | 19 | | 40 017 68 | Bags Brown paper Size: 16 Width: 196mm ; Gusset: 115mm ; Length: 385mm |
| | 20 | | 40 017 70 | Bags Brown paper Size: 20 Width: 210mm ; Gusset: 127mm ; Length: 390mm |
| | 21 | | 40 017 72 | Bags Brown paper Size: 25 Width: 208mm ; Gusset: 127mm ; Length: 430mm |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|--------------------|-----|------|-----------|--|
| PLASTIC BOTTLES | 22 | | 74 110 22 | Plastic Rectangular Medical Bottle with screw cap Capacity: 50ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 23 | | 74 110 41 | Plastic Rectangular Medical Bottle with screw cap Capacity: 100ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 24 | | 74 110 85 | Plastic Rectangular Medical Bottle with screw cap Capacity: 200ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 25 | | 74 110 97 | Plastic Round Medical Bottle with screw cap Capacity: 500ml Colour: Amber Screw cap: 28mm Screw neck: 28mm |
| | 26 | | 74 110 98 | Plastic Round Medical Bottle with Spray Screw Cap Capacity: 500ml Colour: Amber |
| | 27 | | 74 110 99 | Plastic Round Medical Bottle with Pump Screw Cap Capacity: 500ml Colour: Amber |
| | 28 | | 74 130 30 | Plastic Round Bottle with Spray Screw Cap Capacity: 100ml Colour: Opaque White |
| | 29 | | 74 120 21 | Polyethylene bottle with handles and a screw cap Capacity: 2.5l Colour: Opaque White |
| GLASS BOTTLES | 30 | | 74 020 10 | Glass bottle with a dropping pipette screw-top lid Capacity: 10ml |
| | 31 | | 74 050 21 | Glass bottle with a plastic screw-top lid Capacity: 50ml |
| | 32 | | 74 050 31 | Glass bottle with a plastic screw-top lid Capacity: 100ml |
| | 33 | | 74 050 41 | Glass bottle with a plastic screw-top lid Capacity: 200ml |
| | 34 | | 74 050 51 | Glass bottle with a plastic screw-top lid Capacity: 500ml |
| MEDICINE MEASURES | 35 | | 74 020 42 | Squeeze-type dropper bottle with a screw cap Capacity: 15 ml |
| | 36 | | 76 716 11 | Medicine Measure Spoon 5ml Graduated at 2.5ml and 5ml |
| | 37 | | 30 372 01 | Measuring Tumbler Graduated: 5 - 50ml |
| OINTMENT CONTAINER | 38 | | 39 010 19 | Ointment Container with snap-on Cap Capacity: 30g Colour: Clear |
| | 39 | | 39 010 20 | Ointment Container with snap-on Cap Capacity: 50g Colour: White Opaque |
| | 40 | | 39 010 21 | Ointment Container with snap-on Cap Capacity: 100g Colour: White Opaque |
| | 41 | | 74 260 45 | Ointment Container with screw on lid Capacity: 500ml Colour: White Opaque |
| SUNDRIES | 42 | | 30 208 03 | Sputum Cups with lid. ± 150ml |
| | 43 | | 30 092 10 | Denture container ± 250ml |
| | 44 | | 35 025 10 | Urine Specimen collection bottle ± 90ml |
| | 45 | | 30 045 70 | Tongue Depressor Width: ± 17mm; Length: ± 150mm Thickness: ± 1.5mm |
| | 46 | | 30 010 12 | Round Wood Stick Applicator Diameter: 2mm; Length: 150mm |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|-----------|-----|------|-----------|---|
| STORAGE | 47 | | 74 800 21 | Plastic storage bin Size 1 Length: 135mm Width: 105mm Height: 70mm For loads ≤ 3kg |
| | 48 | | 74 800 22 | Plastic storage bin Size 2 Length: 185mm Width: 105mm Height: 70mm For loads ≤ 3kg |
| | 49 | | 74 800 23 | Plastic storage bin Size 3 Length: 275mm Width: 140mm Height: 120mm For loads ≤ 10kg |
| | 50 | | 74 800 24 | Plastic storage bin Size 4 Length: 380mm Width: 205mm Height: 170mm For loads ≤ 15kg |
| | 51 | | 74 800 25 | Plastic storage bin Size 5 Length: 380mm Width: 420mm Height: 170mm For loads ≤ 20kg |
| | 52 | | 39 010 15 | Corrugated Cardboard Carton Length: 305mm Width: 102mm Height: 120mm |
| | 53 | | 39 010 16 | Corrugated Cardboard Carton Length: 305mm Width: 152mm Height: 120mm |
| | 54 | | 39 010 17 | Corrugated Cardboard Carton Length: 300mm Width: 120mm Height: 110mm |
| | 55 | | 39 010 18 | Corrugated Cardboard Carton Length: 390mm Width: 102mm Height: 120mm |

| SURGICAL | | | | |
|--------------------|-----|------|-----------|--|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| SURGICAL SUNDRIES | 1 | | 30 465 10 | Medical Razor |
| | 2 | | 30 248 38 | Stitch Cutter – disposable Length: 9cm - 12cm |
| BLADES AND HANDLES | 3 | | 30 045 25 | Blades Surgical Scalpel Bp No 10 - for making large straight incisions through skin and subcutaneous tissues |
| | 4 | | 30 045 26 | Blades Surgical Scalpel Bp No 11 - elongated triangular blade sharpened along the hypotenuse edge |
| | 5 | | 30 045 27 | Blades Surgical Scalpel Bp No 12 - for suture removal - curved and pointed cutting edge at one end |
| | 6 | | 30 045 28 | Blades Surgical Scalpel Bp No 15 - with a small, curved cutting edge |
| | 7 | | 30 045 29 | Blades Surgical Scalpel Bp No 20 - with a curved cutting edge and a flat, unsharpened back edge |
| | 8 | | 30 045 30 | Blades Surgical Scalpel Bp No 21 - with a large, curved blade |
| | 9 | | 30 045 31 | Blades Surgical Scalpel Bp No 22 - with a curved cutting edge and a flat, unsharpened back edge |
| | 10 | | 30 045 32 | Blades Surgical Scalpel Bp No 23 - "leaf-shaped" blade sharpened along its leading edge. |
| | 11 | | 30 045 33 | Blades Surgical Scalpel Bp No 24 - semi-circular in shape sharpened along its leading edge |
| | 12 | | 60 0039 | Surgical Scalpel Handle No.3 for blade sizes 10 - 17 |
| | 13 | | 60 0040 | Surgical Scalpel Handle No.4 for blade sizes >18 |
| SKIN GRAFT | 14 | | 30 045 15 | Skin Graft Blades: Bodenham Humby 157mm x 18mm x 0.23mm |
| | 15 | | 30 376 03 | Mesh Derma-carrier Size 1.5:1 Length: 200mm Width: 80mm |
| | 16 | | 30 376 02 | Mesh Derma-carrier Size 3:1 Length: 200mm Width: 80mm |
| ENDOSCOPY | 17 | | LapPort1 | Laparoscopic Port Set - 4 Ports: Optical bladeless tip trocars and stability sleeves. 2 Ports: Length: 100mm Diameter: 5mm; 2 Ports: Length: 100mm Diameter: 10mm |
| | 18 | | LapPort2 | Laparoscopic Port - Optical bladeless tip trocar and stability sleeve. Length: 100mm Diameter: 5mm |
| | 19 | | LapPort3 | Laparoscopic Port - Optical bladeless tip trocar and stability sleeve. Length: 100mm Diameter: 10mm |
| STAPLERS | 20 | | LCStap | Linear Cutter Stapler 80mm Width: 3.8mm |
| | 21 | | LCReload | Linear Cutter Stapler Reloads 80mm Width: 3.8mm |
| | 22 | | 30 473 24 | Skin Staples Remover - DISPOSABLE |
| | 23 | | 30 473 25 | Skin Staples Remover - REUSABLE |
| | 24 | | 30 473 26 | Skin Staples Gun – Single use Pre-Packed with 35 staples Width: 5.6mm Height: 3.8mm |

| SURGICAL | | | | |
|--------------|-----|------|-----------|--|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| DIATHERMY | 25 | | 30 249 98 | Diathermy Plates Child D135P |
| | 26 | | 30 249 99 | Diathermy Plates Adult |
| | 27 | | 30 250 14 | Diathermy Hand Control |
| | 28 | | 30 440 30 | Diathermy Scratch Pads Width: 8cm Length: 10cm |
| FEEDING | 29 | | 33 006 69 | Naso-jejunal feeding tube - 6Fr Length: 550mm +/- 10% |
| | 30 | | 33 006 70 | Naso-jejunal feeding tube - 8Fr Length: 550mm +/- 10% |
| | 31 | | 33 006 71 | Naso-jejunal feeding tube – 10Fr Length: 1100mm +/- 10% |
| | 32 | | 33 006 72 | Naso-jejunal feeding tube – 12Fr Length: 1100mm +/- 10% |
| | 33 | | 33 006 79 | Nasogastric tube - Ryles type – 8Fr Length: 100 cm +/- 10% Colour: Blue |
| | 34 | | 33 006 80 | Nasogastric tube - Ryles type – 10Fr Length: 100 cm +/- 10% Colour: Black |
| | 35 | | 33 006 81 | Nasogastric tube - Ryles type – 12Fr Length: 100 cm +/- 10% Colour: White |
| | 36 | | 33 006 82 | Nasogastric tube - Ryles type – 14Fr Length: 100 cm +/- 10% Colour: Green |
| | 37 | | 33 006 83 | Nasogastric tube - Ryles type – 16Fr Length: 100 cm +/- 10% Colour: Orange |
| | 38 | | 33 006 84 | Nasogastric tube - Ryles type – 18Fr Length: 100 cm +/- 10% Colour: Red |
| | 39 | | 33 006 85 | Nasogastric tube - Ryles type – 20Fr Length: 100 cm +/- 10% Colour: Yellow |
| GASTRIC | 40 | | 33 135 17 | Flatus Tube 28Fr Minimum length: 500mm |
| | 41 | | 33 150 13 | Stomach Wash Out Tube 28Fr Minimum Length: 800mm |
| | 42 | | 33 150 17 | Stomach Wash Out Tube 30Fr Minimum Length: 800mm |
| | 43 | | 33 150 19 | Stomach Wash Out Tube 32Fr Minimum Length: 800mm |
| | 44 | | Sengs | Sengstaken Tube 21Fr Length: 100cm |
| WOUND DRAINS | 45 | | 30 226 10 | Corrugated Drains - 350mm x 45mm (±5%) |
| | 46 | | 30 226 11 | Corrugated Drains - 350mm x 110mm (±5%) |
| | 47 | | 30 226 12 | Corrugated Drain Sheet with drainage lumens Size: 40cm x 110mm (±5%) |
| | 48 | | PD1 | Pencil Drains - 6mm x 300mm |
| | 49 | | PD2 | Pencil Drains - 10mm x 300mm |
| | 50 | | PD3 | Pencil Drains - 20mm x 300mm |
| | 51 | | WDR1 | Wound Drain - Hubless 10 FR Round, 1/8" Trocar |
| | 52 | | WDR2 | Wound Drain - Hubless 15 FR Round, 3/16" Trocar |
| | 53 | | WDR3 | Wound Drain - Hubless 19 FR Round, 1/4" Trocar |
| | 54 | | WDR4 | Wound Drain - Hubless 19 FR Round, 1/4" Bendable Trocar |
| | 55 | | WDR9 | Wound Drain - Hubless, without Trocar 24Fr (8.0mm) Round |
| | 56 | | WDF2 | Wound Drain - 10mm Flat, 3/4 Fluted, 3/16" Trocar |
| | 57 | | WDF3 | Wound Drain - 10mm Flat, Full Fluted |

SURGICAL

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|-------------------|------------|-------------|--------------|---|
| | 58 | | WDF4 | Wound Drain - 10mm Flat, Full Fluted, 3/16" Trocar |
| RESERVOIRS | 59 | | CWDR1 | Closed wound drainage reservoir 150ml |
| | 60 | | CWDR2 | Closed wound drainage reservoir 300ml |
| | 61 | | CWDR3 | Closed wound drainage reservoir 450 - 500ml |
| | 62 | | CWDR4 | Bulb Suction Reservoir 100 ml |

| OPHTHALMOLOGY | | | | |
|----------------------|-----|------|-----------|---|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| SHIELDS | 1 | | ESheildA | Eye Shield - Adult 6.2cm x 7.5cm |
| | 2 | | ESheildP | Eye Shield - Paediatric 5.1cm x 4.1cm |
| CONFORMERS | 3 | | EConf1 | Eye Conformer Size: 19mm x 22mm |
| | 4 | | EConf2 | Eye Conformer Size: 21mm x 24mm |
| | 5 | | EConf3 | Eye Conformer Size: 22mm x 25mm |
| ORBITAL BALL IMPLANT | 6 | | Orb1 | Orbital Ball Implant Size: 10-14mm |
| | 7 | | Orb2 | Orbital Ball Implant Size: 16mm |
| | 8 | | Orb3 | Orbital Ball Implant Size: 18mm |
| | 9 | | Orb4 | Orbital Ball Implant Size: 20mm |
| OPHTHALMIC KNIVES | 10 | | CresKnife | Ophthalmic Crescent Knife |
| | 11 | | BlebKnife | Ophthalmic Bleb Knife |
| | 12 | | SlitKnife | Ophthalmic Slit Knife |
| | 13 | | StabKnife | Ophthalmic Stab Knife |
| OPHTHALMIC | 14 | | Spear | Ophthalmic Spear Including Stick |
| | 15 | | Trypan | Trypan Blue Ophthalmic Solution |
| | 16 | | MonStent | Mini Monoka Stent Length: 40mm ID: 0.3mm OD: 0.64mm |
| | 17 | | AValve | Glaucoma Drainage Device - Ahmed Valve |

| CSSD | | | | |
|------------|-----|------|-----------|--|
| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
| TAPE | 1 | | 30 577 49 | Autoclave Indicator Tape Size:12mm x 50m |
| | 2 | | 30 577 50 | Autoclave Indicator Tape Size: 24mm x 50m |
| | 3 | | 30 578 05 | Autoclave Adhesive Tape Size: 24mm x 50m |
| | 4 | | CCTWhite | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: White |
| | 5 | | CCTBlue | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Blue |
| | 6 | | CCTGreen | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Green |
| | 7 | | CCTRed | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Red |
| | 8 | | CCTYellow | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Yellow |
| | 9 | | CCTPurple | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Purple |
| | 10 | | CCTOrange | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Orange |
| | 11 | | CCTBrown | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Brown |
| | 12 | | CCTBlack | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Black |
| INDICATORS | 13 | | 30 338 35 | Steam Sterilisation Indicator Strips Size: 1.5cm x 20cm |
| | 14 | | 30 338 32 | Bowie Dick Indicator for Autoclaves |
| | 15 | | 30 338 30 | Ethylene Oxide (EO) Chemical Indicator Strip (Attest strip) Size: 1.5cm x 20cm |
| | 16 | | BioTT | Biological Indicator Test Tube (BT220 equivalent) |
| | 17 | | BioPP | Biological Indicator Thermal Printing Paper |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|---------------------|-----|------|-----------|---|
| STERILISATION BAGS | 18 | | 30 018 10 | Paper Bag for Steam Sterilisation Size: 1 Width: 90mm Length: 100mm Gusset: 50mm |
| | 19 | | 30 018 11 | Paper Bag for Steam Sterilisation Size: 2 Width: 90mm Length: 230mm Gusset: 50mm |
| | 20 | | 30 018 12 | Paper Bag for Steam Sterilisation Size: 3 Width: 115mm Length: 190mm Gusset: 35mm |
| | 21 | | 30 018 13 | Paper Bag for Steam Sterilisation Size: 4 Width: 125mm Length: 260mm Gusset: 50mm |
| | 22 | | 30 018 15 | Paper Bag for Steam Sterilisation Size: 6 Width: 190mm Length: 320mm Gusset: 60mm |
| | 23 | | 30 018 16 | Paper Bag for Steam Sterilisation Size: 7 Width: 180mm Length: 380mm Gusset: 75mm |
| | 24 | | 30 018 17 | Paper Bag for Steam Sterilisation Size: 8 Width: 65mm Length: 480mm Gusset: 30mm |
| | 25 | | 30 018 18 | Paper Bag for Steam Sterilisation Size: 9 Width: 120mm Length: 135mm Gusset: 30mm |
| | 26 | | 30 018 19 | Paper Bag for Steam Sterilisation Size: 11 Width: 290mm Length: 445mm Gusset: 125mm |
| | 27 | | 30 018 20 | Paper Bag for Steam Sterilisation Size: 12 Width: 330mm Length: 560mm Gusset: 150mm |
| | 28 | | 30 018 31 | Paper Bag for Steam Sterilisation Size: 1020 Width: 65mm Length: 200mm Gusset: 35mm |
| STERILISATION PAPER | 29 | | 30 445 32 | Crêped Sterilisation Paper White/Blue/Green Size: 300mm x 300mm |
| | 30 | | 30 445 31 | Crêped Sterilisation Paper White/Blue/Green Size: 600mm x 600mm |
| | 31 | | 30 445 33 | Crêped Sterilisation Paper White/Blue/Green Size: 600mm X 900mm |
| | 32 | | 30 445 36 | Crêped Sterilisation Paper White/Blue/Green Size: 900mm x 900mm |
| | 33 | | 30 445 35 | Crêped Sterilisation Paper White/Blue/Green Size: 900mm x 1300mm |
| | 34 | | 30 445 34 | Crêped Sterilisation Paper White/ Blue/Green Size: 900mm x 1500mm |
| | 35 | | 30 445 29 | Crêped Sterilisation Paper White/ Blue/Green Size: 1000mm x 1400mm |
| | 36 | | 30 445 27 | Crêped Sterilisation Paper White/Blue/Green Size: 1500mm x 1500mm |
| STERILISATION WRAP | 37 | | 39 127 49 | Sterilisation Wrap (Spunbound) - Regular Size: 600mm x 600mm |
| | 38 | | 39 127 50 | Sterilisation Wrap (Spunbound) - Regular Size: 900mm x 900mm |
| | 39 | | 39 127 51 | Sterilisation Wrap (Spunbound) - Regular Size: 1200mm x 1200mm |
| | 40 | | 39 127 52 | Sterilisation Wrap (Spunbound) - Regular Size: 1400mm x 1400mm |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: |
|-----------------------|-----|------|-----------|--|
| STERILISATION REELS | 41 | | 30 445 38 | Flat Sterilising Reels without gusset Size: 50mm x 200m |
| | 42 | | 30 445 39 | Flat Sterilising Reels without gusset Size: 75mm x 200m |
| | 43 | | 30 445 40 | Flat Sterilising Reels without gusset Size: 100mm x 200m |
| | 44 | | 30 445 41 | Flat Sterilising Reels without gusset Size: 150mm x 200m |
| | 45 | | 30 445 42 | Flat Sterilising Reels without gusset Size: 200mm x 200m |
| | 46 | | 30 445 43 | Flat Sterilising Reels without gusset Size: 250mm x 200m |
| | 47 | | 30 445 45 | Flat Sterilising Reels with gusset Size: 100mm x 150mm |
| | 48 | | 30 445 46 | Flat Sterilising Reels with gusset Size: 100mm x 200mm |
| | 49 | | 30 445 47 | Flat Sterilising Reels with gusset Size: 100mm x 250mm |
| DETERGENTS & CLEANERS | 50 | | 40 030 57 | Detergent Liquid Cleaner for Theatre instruments (Endozime equivalent) |
| | 51 | | 30 073 10 | Enzymatic Presoak & Cleaner - 3.8litres (Prolystica Equivalent) |
| | 52 | | 30 073 11 | Neutral Detergent - 3.8litres |
| | 53 | | 30 073 12 | Alkaline Detergent- 3.8litres |
| | 54 | | 60 00 09 | Disinfectant Sterilising Solution (Cidex equivalent) |
| | 55 | | SurgiSt | Surgical Instrument rust and scale remover (Surgistain equivalent) |
| | 56 | | NBrush | Nylon Instrument Brush |
| | 57 | | SSBrush | Stainless Steel Instrument Brush |
| CSSD PROTECTION | 58 | | EyeWash1 | Sterile Eye Wash Solution 200- 250ml |
| | 59 | | EyeWash2 | Sterile Eye Wash Solution 500ml |
| | 60 | | HRApron | Heat-resistant Apron |
| | 61 | | HRGloves | Heat-resistant Autoclave Gloves |
| | 62 | | NitGlove1 | Non-Sterile Nitrile Decontamination Gloves Size: Small |
| | 63 | | NitGlove2 | Non-Sterile Nitrile Decontamination Gloves Size: Medium |
| | 64 | | NitGlove3 | Non-Sterile Nitrile Decontamination Gloves Size: Large |
| | 65 | | NitGlove4 | Non-Sterile Nitrile Decontamination Gloves Size: X-Large |

PERSONAL IDENTIFICATION BANDS

| PERSONAL IDENTIFICATION BAND - ADULT | | |
|--|--|--|
| Manufactured from soft, pliable, durable, and hypoallergenic plastic or other suitable material. Band width 2cm to allow for printing, in an indelible manner, the patient's name and other particulars, Band to be supplied with a locking device which cannot ordinarily be opened after first closure. Locking device from 12cm - 25cm Packaging: Box of 100; Unit = 1 Box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|----------------|--------------|---|
| | 30 028 07 | Personal Identification Band - Adult Yellow Locking device from 12cm - 25cm |
| | 30 028 08 | Personal Identification Band - Adult Pink Locking device from 12cm - 25cm |
| | 30 028 09 | Personal Identification Band - Adult White Locking device from 12cm - 25cm |

| PERSONAL IDENTIFICATION BAND - CHILD | | |
|---|--|--|
| Manufactured from soft, pliable, durable, and hypoallergenic plastic or other suitable material. Band width 2cm to allow for printing, in an indelible manner, the patient's name and other particulars. Band to be supplied with a locking device which cannot ordinarily be opened after first closure. Locking device from 9cm - 21cm Packaging: Box of 100; Unit = 1 Box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|----------------|--------------|--|
| | 30 028 10 | Personal Identification Band - Child White Locking device from 9cm - 21cm |
| | 30 028 12 | Personal Identification Band - Child Pink Locking device from 9cm - 21cm |
| | 30 028 13 | Personal Identification Band - Child Yellow Locking device from 9cm - 21cm |

| PERSONAL IDENTIFICATION BAND - NEONATES | | |
|--|--|--|
| Manufactured from extra-soft , pliable, durable, and hypoallergenic plastic or other suitable material – to prevent abrasion of delicate skin. Band width 2cm to allow for printing, in an indelible manner, the patient's name and other particulars. Band to be supplied with a locking device which cannot ordinarily be opened after first closure. Plastic press studs must fasten with ease. Locking device from 6cm - 12cm Packaging: Box of 100; Unit = 1 Box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|----------------|--------------|--|
| | 30 028 11 | Personal Identification Band - Neonate and Premature Babies White Locking device from 6cm - 12cm |

DISPOSABLE UNDERWEAR

| DISPOSABLE UNDERWEAR | | |
|---|-----------|--|
| Manufactured from non-transparent non-woven elasticised yarn ≥ 35 g/m ² . The leg and waist opening must be over-locked and must have an elasticated waist and leg. Packaging: Bale of 100; Unit = 1 bale of 100 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 39 054 13 | Disposable underwear - Small Briefs Waist circumference: 80-90cm Leg circumference: 28-60cm Waist-crotch: ≥ 30 cm |
| | 39 054 14 | Disposable underwear - Medium Briefs Waist circumference: 90-100cm Leg circumference: 32-64cm Waist-crotch: ≥ 32 cm |
| | 39 054 15 | Disposable underwear - Large Briefs Waist circumference: 100-120cm Leg circumference: 36-72cm Waist-crotch: ≥ 35 cm |
| | 39 054 16 | Disposable underwear - X-Large Briefs Waist circumference: 120-140cm Leg circumference: 36-76cm Waist-crotch: ≥ 40 cm |
| | 39 054 17 | Disposable underwear - XX-Large Briefs Waist circumference: 130cm-150cm Leg circumference: 40cm-80cm Waist-crotch: 44cm |
| | 39 054 18 | Disposable underwear - XXX-Large Briefs Waist circumference: 140cm-160cm Leg circumference: 44cm-88cm Waist-crotch: 46cm |

FLUSHABLE COVERS

| FLUSHABLE COVERS | | |
|---|-----------|--|
| Flushable covers made from Toilet Paper 24G/Square Metre Packaging: Bale of 100; Unit = 1 bale of 100 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 39 050 11 | Bedpan Cover - Flushable Size: 30.5cm x 7.6cm x 40.6cm |
| | 39 050 12 | Urinal Covers - Flushable Size: 15.22cm x 7.6cm x 30.5cm |

PROTECTIVE COVERS

| ULTRASOUND PROBE COVER | | |
|---|----------|--|
| Purpose: Sterile plastic sleeve for ultrasound use To be made from ultra-thin polyurethane film. Must not interfere with the view or create artefacts. Must fit both linear and curvilinear probes and the kit must include sterile ultrasound gel and elastic bands. Latex free and sterile Packaging: Pack of 50 Covers | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | USPcover | Ultrasound probe cover Size: ± 10 cm x ± 60 cm |

| TEMPERATURE PROBE COVER | | |
|--|-----------|--|
| Slim fitting condom to cover internal temperature probes. Must fit both adult and paediatric probes. Must allow for accurate temperature measurement without insulation. Must be latex free and hypoallergenic Packaging: Pack of 50 Covers | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 461 10 | Temperature probe cover Length: 15-20cm |

WIPES

| ALCOHOL SWABS | | |
|---|-----------|------------------------------------|
| Purpose: For skin preparation prior to injection or IVI cannulation Composed of viscose rayon non-woven fabric impregnated with minimum of 0,5ml Isopropyl Alcohol 70% in water. Minimum dimensions: 24mm x 30mm Individually packed in aluminium foil pouches. Sterile Packaging: Box of 200 swabs | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 568 79 | Alcohol Swabs - 24mm x 30mm |

| WIPES TO REMOVE ADHESIVENESS FROM THE SKIN | | |
|---|---------|---|
| Purpose: Wipes to remove adhesiveness residue from the skin Individually wrapped sterile wipe Size: 40cm x 40cm Non-irritating and non-sensitizing wipe that moisturises and conditions the skin without greasiness. Must gently clean hydrocolloid, acrylic and rubber-based adhesive residues from the skin without irritation, discomfort, or trauma. Safe for ostomies and geriatric patients. Single use Packaging: Box of 50 wipes | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | AdhWipe | Wipe to remove adhesiveness from the skin. Size: 40cm x 40cm |

GELS

| ECG/ ULTRASOUND/ DEFIBRILLATOR GEL | | |
|---|-----------|--|
| Purpose: Electrically conductive gel composition for use in establishing a low resistance contact between an electrode and skin Medical quality clear, aqueous GEL with the CE mark Must provide extremely high conductive effect and clear, accurate recording without electrical artefacts. Must be highly viscous and resistant to drying. Must be chemically neutral with no interactions with the patient skin or electrodes. Bottle lids must be secure and not pop off when the bottle is squeezed. Packaging: 10 bottles / 10 jerry cans | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 344 30 | Ultrasound Gel with dispenser 200ml |
| | 30 344 32 | Ultrasound Gel Jerry Can 5 Litre |

| LUBRICATING JELLY (KY JELLY) | | |
|---|--------|--|
| Purpose: Jelly for lubricating during gynaecological procedure; electrotherapy; cystoscopy and general lubrications. Non Petroleum based lubricant must be suitable for use on all silicone and PVC tubes. Must be water soluble, Non-greasy, Non-sticky and clear. Sterile sachets Packaging: 50 Sachets | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | LubGel | Lubricating (KY) Jelly 2.5g Sachets |

| ANAESTHETIC GEL | | |
|---|-----------|-------------------------------------|
| Purpose: To anaesthetise the airway or urogenital tract Lubricating jelly with 2g lignocaine hydrochloride per 100g Packaging: Per Unit | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 53 026 10 | Anaesthetic Gel 2% 20ml Tube |

ECG ELECTRODES

SOLID GEL ECG ELECTRODE

Purpose: For diagnostic and monitoring purposes of an ECG tracing for up to 24 hours.

Consists of a **pre-gelled Silver/Silver-Chloride solid gel, spongeless electrode** with a press stud that fits standard ECG leads.

The electrode must be easy to strip from the backing liner.

The adhesive backing must be **hypoallergenic** and **latex free** and provide good contact and resist electrode contact loss due to perspiration.

There must be a non-adhesive edge on one side of the backing to facilitate removal off the patient.

The electrode must be disposable using standard waste disposal procedures.

It must have a shelf-life of 2 years if unopened and be SABS 1461 or CE approved.

Packaging to clearly state: Type and size of electrode ; SABS/CE compliance; Expiry date

Packaging: 50/Package

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 248 03 | Solid Gel ECG Electrode - Adult |
| | 30 248 04 | Solid Gel ECG Electrode - Paediatric < 50mm wide |

SPONGE ECG ELECTRODES

Purpose: For diagnostic and monitoring purposes of an ECG tracing for up to 24 hours.

Consists of a **pre-gelled Silver/Silver-Chloride self-adhering sponge** electrode with a press stud that fits standard ECG leads.

The electrode must be easy to strip from the backing liner.

The sponge backing materials must protect the sensor and gel from surgical and cleaning fluids and still permit maximum flexibility.

There must be a non-adhesive edge on one side of the backing to facilitate removal off the patient.

The electrode must be disposable using standard waste disposal procedures.

It must have a shelf-life of 2 years if unopened and be SABS 1461 or CE approved.

Packaging to clearly state: Type and size of electrode ; SABS/CE compliance; Expiry date

Packaging: 50/Package

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 248 09 | Sponge ECG Electrode - Paediatric < 50mm wide |
| | 30 248 18 | Sponge ECG Electrode - Adult |

TAPE/CLOTH ECG ELECTRODE

Purpose: For diagnostic and monitoring purpose of an ECG tracing > 24 hours

Consists of a **pre-gelled Silver/Silver-Chloride self-adhering tape/cloth** electrode with a press stud that fits standard ECG leads.

The electrode must be easy to strip from the backing liner.

The sponge backing materials must protect the sensor and gel from surgical and cleaning fluids and still permit maximum flexibility.

There must be a non-adhesive edge on one side of the backing to facilitate removal off the patient.

The electrode must be disposable using standard waste disposal procedures.

It must have a shelf-life of 2 years if unopened and be SABS 1461 or CE approved.

Packaging to clearly state: Type and size of electrode ; SABS/CE compliance; Expiry date

Packaging: 50/Package

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 248 27 | Tape ECG Electrode - Paediatric < 50mm wide |
| | 30 248 37 | Tape ECG Electrode - Neonate < 30mm wide |
| | 30 248 38 | Tape ECG Electrode - Adult |

PLASTIC APRONS

| | | |
|--|--|--|
| FULL BODY APRON | | |
| FULL BODY APRON | | |
| <p>Purpose: To be protect clothing from bodily fluids or other organic product A single use plastic full body apron Must be noiseless and made from smooth waterproof plastic material. Ties length: ± 50cm; Ties Width: ± 5cm; Width: ± 65cm; Thickness: 25 micron</p> <p>For single use only Packaging: Bale of 50</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 39 010 15 | Full Body Apron Medium - General use Length from neck: ± 110cm Colour: Blue |
| | 39 010 16 | Full Body Apron Large - General use Length from neck: ± 130cm Colour: Blue |
| | 39 010 18 | Full Body Apron Medium - Sluice room Length from neck: ± 110cm Colour: Red |
| | 39 010 19 | Full Body Apron Large - Sluice room Length from neck: ± 130cm Colour: Red |
| | 39 010 21 | Full Body Apron Medium - Infectious areas Length from neck: ± 110cm Colour: Yellow |
| | 39 010 22 | Full Body Apron Large - Infectious areas Length from neck: ± 130cm Colour: Yellow |
| | 39 010 24 | Full Body Apron Medium - Clean areas (kitchen) Length from neck: ± 110cm Colour: White |
| | 39 010 25 | Full Body Apron Large - Clean areas (kitchen) Length from neck: ± 130cm Colour: White |

| | | |
|--|--|--|
| FULL BODY APRON - HEAVY DUTY | | |
| <p>Purpose: To be protect clothing from bodily fluids or other organic product</p> <p>Made from soft reusable semi-transparent, waterproof polyvinyl chloride PVC Gauge: 6 - 8 with waterproof seams. Must have soft, adjustable neck bands and tapes for waist. Reusable Packaging: Per Apron</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 39 010 26 | Surgical Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Light Green |
| | 54 020 10 | Laundry Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Red |

NON-DISPOSABLES

STETHOSCOPE - HIGH END

Purpose: To identify, listen to, and study heart, lung, and other body sounds for physical assessment and diagnosis
Consisting of a dual-sided chest piece that provides high acoustic sensitivity for both high and low frequency sounds with both a floating diaphragm and open bell.

Tubing must be flexible, resilient and phthalate plasticizers free, **± 70cm** long.

Headset tension must be adjustable for various head sizes.

The ear tips must provide a good acoustic seal and comfortable fit and the ear tubes must be angled to align with the ear canals and must be able to rotate.

All components must be latex free.

A **5-year warranty** must be included. **Per stethoscope**

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 542 06 | Stethoscope – Infant Diaphragm: ± 2.7cm |
| | 30 542 07 | Stethoscope – Paediatric Diaphragm: ± 3.3cm |
| | 30 542 08 | Stethoscope – Adult Diaphragm: ± 5cm |

DIGITAL THERMOMETER

Purpose: Single patient use thermometers to record oral and axillary temperatures

Must be made of BPA free plastic with a flexible tip. Must have a storage case.

Must be capable of giving fast accurate reading within 10 seconds.

Must cover temperature ranges of at least: **34 °C to 42°C**. (Not Fahrenheit)

Must have a memory to record the last temperature reading.

Must have automatic shut off function; Must be supplied with a battery.

For oral and axillary use.

Must be completely sealed and compatible with alcohol disinfectant.

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 244 00 | Thermometer – Digital temperature ranges of at least: 34°C to 42°C . |

PRESSURE BAG - NON-DISPOSABLE

Purpose: To apply pressure on infusion bags and speed up intravenous fluid administration

Made from transparent high-grade medical plastic to allow a clear view of the fluid bag and the fluid level.

Must have a wrap-around design and provide consistent pressure on the complete fluid bag.

The wrap-around design must be easy to open with a hook for hanging the fluid bag.

Replacing the fluid bag must be quick and easy.

The pressure bag must have a shock resistant manometer that indicates precise pressure with a range from **0-300mmHg**.

The hand inflator must allow large volumes and quick inflation.

A **3-way stop cock** must allow for quick deflation.

The pressure bag must be reusable, durable, easy to clean and latex free.

Pressure bag to be guaranteed to last 5 years.

Packaging: each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--------------------------------------|
| | ND-PB | Pressure Bag – Non-Disposable |

MANOMETER – ENDOTRACHEAL AND SUPRAGLOTTIC CUFF PRESSURE MONITOR

Purpose: to inflate and monitor the pressure of laryngeal masks and intubation tubes with low-pressure cuffs to reduce the risk of pressure necrosis and mucosal ischemia as well as preventing the risk of aspiration.

Consists of a large dial with easy-to-read display graded in cmH₂O.

The display must clearly indicate the ideal ranges for either tracheal tubes or laryngeal tubes.

22–32 cmH₂O for tracheal tubes

32–60 cmH₂O For Laryngeal Tubes and Laryngeal Masks

A luer-connection must allow the gauge to be connected to laryngeal masks or intubation tubes, either directly or via the connecting tube.

Must have a release valve to allow air release in the case of over-inflation.

The cuff must be inflated using a sensitive and easy to compress inflation bulb.

Must have a suspension hook for easy rail connection.

The manometer must be ergonomic, durable, and easy to clean.

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | ET-M | Manometer - Endotracheal and Supraglottic Cuff Pressure Monitor |

PASTEURISATION

| | | |
|---|--|--|
| PASTEURISATION BOTTLE | | |
| <p>Purpose: Reusable glass bottles suitable for the pasteurisation of human donor milk</p> <p>Made from glass that is heat-resistant and autoclavable up to 120°C. Volume: 125ml Neck closure: 48mm single thread Must be supplied with a suitable lid made of metal with plastisol that is also autoclavable. Unit = 1 bottle</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 047 00 | Pasteurisation Bottle - Glass. Volume: 125ml |

| | | |
|--|--|--|
| LID FOR GLASS PASTEURISATION BOTTLE | | |
| <p>Purpose: Lid to seal pasteurisation jars and prevent air ingress or leakage Metal Plastisol Lid with 48mm single thread Heat-resistant and autoclavable up to 120°C Unit = 1 lid</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 047 01 | Lid for Glass Pasteurisation Bottle 48mm single thread |

| | | |
|--|--|--|
| LID FOR GLASS PASTEURISATION BOTTLE | | |
| <p>Purpose: Lid to seal pasteurisation jars and prevent air ingress or leakage Metal Plastisol Lid with 48mm single thread Heat-resistant and autoclavable up to 120°C Unit = 1 lid</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 047 01 | Lid for Glass Pasteurisation Bottle 48mm single thread |

| | | |
|--|--|--|
| LIDS WITH FOIL SEAL FOR HUMAN MILK STORAGE BOTTLES - REUSABLE | | |
| <p>Purpose: Reusable screw cap lids with foil seal to provide secure sealing of BPA-free Pasteurisation Bottles</p> <p>Consists of a plastic screw cap that fits securely on the BPA-free plastic bottle. The screw cap must have a foil seal to provide secure sealing during pasteurisation, and transportation up to the point of use. No water should leak into the human milk during pasteurisation and cooling. The seal should have a tab to avoid touching the neck of the bottle when removing the tamper evident foil seal. Heat resistant and autoclavable up to 120 °C. Packaging: Each</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 047 03 | Lids with foil seal for Human Milk Storage Bottles - Reusable |

| | | |
|---|--|--|
| BOTTLE BRUSH: | | |
| <p>Purpose: For the cleaning of pasteurisation bottles Length including brush: 40cm. Brush: 10cm, Handle: 30cm Flexible bristles: 5cm. Handle to have a loop for hanging. Packaging: Each</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 052 09 | Bottle Brush: Length including brush: 40cm. Brush: 10cm, Handle: 30cm Flexible bristles: 5cm |

| | | |
|---|--|--|
| SMALL-SCALE PASTEURISER LABELS | | |
| <p>Purpose: Small-scale pasteuriser labels that fit a Dymo Label writer for the labelling of pasteuriser containers</p> <p>The labels must be water and moisture safe. They must resist degradation in cold storage and high temperatures. They must fit a Dymo label writer and come as a roll with 130 labels. Label size: 28mm x 89mm</p> | | |

Unit = 2 Rolls of 130 labels

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|---|
| | PL1 | Small-scale Pasteuriser Labels Size: 28mm x 89mm |

PLASTICISED PASTEURISER LABELS

Purpose: Medium/Large-scale pasteuriser labels for the labelling of pasteuriser bottles

Made from waterproof **Polyolefin** that is easy to peel off and on plastic or glass bottles.
The labels must remain intact even with immersion in 62° water.
Labels must be blank and allow for use with a permanent marker. Must be easy to remove from the sheet.

Unit = 1 Roll of 2440 labels

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | PL2 | Plasticised Pasteuriser Labels Size: 100mm x 25mm |

NEONATAL

LITMUS PAPER

Litmus Paper to indicate alkali or acidic medium

Packaging: Booklet with 20 strips; Unit = 1 booklet

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|-------------------------------|
| | 30 445 23 | Litmus Paper Dark Blue |
| | 30 445 24 | Litmus Paper Red |

pH INDICATOR TEST STRIPS

Purpose: To indicate if NGT's are in the stomach

Sensitive pH test strips with a wide range of readings from pH **4.5-9.0**

Packaging: Each; Unit = 1 strip

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---------------------------------|
| | 30 445 25 | pH Indicator Test strips |

PHOTOTHERAPY EYE-SHIELD

Purpose: To protect patients' eyes during phototherapy while allowing phototherapy light to reach the baby's head to give the maximum treatment effect.

Consists of a wraparound one-piece head strap with soft, padded eye shields that fit over the eyes.

The head strap must be secured with adjustable velcro straps that must be independently adjustable.

Must fit securely yet be easy to move or adjust when necessary.

Material must be soft, hypoallergenic, latex and DEHP free.

Must offer sufficient ventilation to the eye and not compress the eye or nasal bridge when secured.

Individually packed

Packaging: Unit = 10 Shields

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 460 51 | Phototherapy Eye-Shield - Preemie, head circumference \pm 26-32cm |
| | 30 460 52 | Phototherapy Eye-Shield - Regular, head circumference \pm 33-38 cm |

UMBILICAL CORD CLAMP

Clamp must be constructed for easy application in wet, slippery conditions.

Serrations and stoppers must prevent the clamp from slipping but must not cut through the cord. Clamp must remain securely locked after application.

Individual sterile pack. Gamma – sterilised. CKS 312

Packaging: Packet of 100; Unit = 1 packet of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--------------------------------------|
| | 30 072 01 | Umbilical Cord Clamp. \pm 5cm long |

AMNIOTOMY HOOK

Purpose: As a tool for the artificial rupture of membranes - to induce and impulse the process of delivery.

Must perforate the membranes without trauma to the surrounding tissue.

Length: **25 -30cm** Width: **1-2cm** tapering off to **0.5cm** nearing the hook; Curve on hook: **0.3mm**.

Curved to ensure a positive grip in the user's hand and to follow the shape of the pelvic outlet.

Latex free, Sterile, and individually packed in a peel pouch that is easy to open.

For single use only

The following must be noted on the packaging:

Trade name; Manufacturing site; Expiry date; CE accredited.

Packaging: Unit = 1 packet of 50

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|------------------------------|
| | 30 326 10 | Amniotomy hook - 25cm |

CONNECTORS

| TUBING CONNECTORS | | |
|--|--|--|
| <p>Purpose: For connecting different sizes of gas / fluid tubing Tough, shatterproof, plastic connectors with barbed tapered ends to help prevent disconnection. Manufactured from: Vinyl, Polypropylene, or Polyethylene Latex free, Clinically clean, autoclavable Packaging: Packet of 50</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 089 23 | Y-Shaped Tubing Connector - Large OD: ±10mm All sides must be of equal length – 30mm length from point to junction |
| | 30 089 36 | Straight Tubing Connector - Open-ended . Double tapered Small OD: 7-11mm; ID: 5mm; Total length: 6cm |
| | 30 089 37 | Straight Tubing Connector - Open-ended . Double tapered Large OD: 10 - 14mm ID: 8mm Total length: 6cm |

| SPIGOTS | | |
|---|--|--|
| <p>Purpose: Solid plastic, tough, shatterproof, tapered from specified diameter to point. Tough, shatterproof, plastic connectors Manufactured from: Vinyl, Polypropylene, or Polyethylene Latex free, Clinically clean, autoclavable Packaging: Packet of 50</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 534 80 | Spigot - Small Diameter: 9mm ; Length: ±40mm |
| | 30 534 82 | Spigot - Large Diameter: 10mm ; Length: ± 52mm |

DRESSING PACKS

| DRESSING PACK FOR MINOR SURGICAL PROCEDURES | | |
|--|--|--|
| <p>Purpose: Sterile Surgical tray - pre-packed for minor surgical procedures Must contain 1 x tray $\pm 175\text{mm} \times 150\text{mm}$ with 2 compartments. 1st compartment $\pm 110\text{mm} \times 40\text{mm}$, depth $\pm 40\text{mm}$; 2nd compartment $\pm 110\text{mm} \times 100\text{mm}$, depth $\pm 40\text{mm}$ The tray must have a $\pm 10\text{mm}$ width seal and must open with relative ease. The tray must have rounded corners with no serrated edges.</p> <p>Packaging to clearly state: Tray size and components; Expiry date Packaging: Box of 10 trays; Unit = 1 box of 10</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 432 15 | <p>Large Tray Dressing Pack for Minor Surgical Procedures The pack must contain the following: 1 x Plastic disposable bag $\pm 400\text{mm} \times 250\text{mm}$; 2 x Protection sheet $\pm 400\text{mm} \times 400\text{mm}$ 2 x Medical towel, serviette $\pm 330\text{mm} \times 300\text{mm}$ 5 x Sterile 8 ply Gauze swabs plain $\pm 100\text{mm} \times 100\text{mm}$; 5 x Sterile Cotton Balls grade 1 2 x Latex gloves, medium 1 x Plastic non-toothed forceps $\pm 130\text{mm}$</p> |
| | 30 432 16 | <p>Multi-Pack Tray Dressing Pack for Minor Surgical Procedures The pack must contain the following: 1 x Plastic disposable bag $\pm 400\text{mm} \times 250\text{mm}$ 2 x Medical towel, serviette $\pm 330\text{mm} \times 300\text{mm}$; 1 x Water drape $\pm 400\text{mm} \times 400\text{mm}$ 2 x Latex gloves, medium 5 x Sterile 8 ply Gauze swabs plain $\pm 100\text{mm} \times 100\text{mm}$; 5 x Sterile Cotton Balls grade 1</p> |
| | 30 432 17 | <p>Catheter Tray Dressing Pack The pack must contain the following: 1 x Plastic disposable bag $\pm 400\text{mm} \times 250\text{mm}$ 2 x Medical towel, serviette $\pm 330\text{mm} \times 300\text{mm}$ 2 pair x Powder free gloves, medium 1 x Crêpe sheet $\pm 500 \times 480\text{mm}$ and Aperture $\pm 150\text{mm} \times 80\text{mm}$ 1 x Fenestrated drape $\pm 500\text{mm} \times 450\text{mm}$ 2 x Protection sheet $\pm 400\text{mm} \times 400\text{mm}$ 1 x 10 ml syringe; 1 x 21G needle ; 1 x 10 ml sterile water 10 x Sterile Cotton Balls grade 1</p> |
| | 30 432 18 | <p>Maternity Swabbing Tray The pack must contain the following: 1 x Plastic disposable bag $\pm 400\text{mm} \times 250\text{mm}$ 2 x Medical towel, serviette $\pm 330\text{mm} \times 300\text{mm}$ 1 pair x Powder free gloves, medium 2 x Protection sheet $\pm 400\text{mm} \times 400\text{mm}$ 1 x Maternity, sanitary towel 10 x Cotton Balls grade 2</p> |

COTTON WOOL AND CLEANSING SWABS

| | | |
|---|--|--|
| COTTON WOOL ROLL | | |
| Purpose: To disperse and mop up after an injection | | |
| Absorbent cotton wool, non-interleaved To comply with the latest issue of SANS 149 | | |
| Packaging: Roll | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 645 05 | Cotton wool roll – 1st Grade 50g roll |
| | 30 645 06 | Cotton wool roll – 2nd Grade 50g roll |

| | | |
|--|--|--|
| COTTON WOOL ABSORBENT BALLS | | |
| Purpose: To disperse and mop up after an injection | | |
| 100% Pure cotton wool, naturally soft and absorbent White, undyed, grade 1; Autoclavable To comply with the latest issue of SANS 149 | | |
| The following must be noted on the packaging and box: Trade name; Size and specification | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 568 03 | Cotton Wool Absorbent Balls – Sterile 5 balls in an easy to open peel pouch Packaging: 1 unit = 100 Pouches |
| | 30 568 04 | Cotton Wool Absorbent Balls – Unsterile Box of 1000 balls |

| | | |
|---|--|--|
| CLEANSING BALL SWABS | | |
| Purpose: To clean surgical or procedural sites | | |
| Made from 1 st grade absorbent cotton wool in a single layer of absorbent gauze. The covering gauze must be secured with cotton thread to prevent the cotton wool from spilling out of it. The swabs must be lightly bound together in bundles of five by means of an easy to break cotton thread. Radio-opaque line to be embedded in the gauze covering. Ball bundles in individual pouches in a labelled box To comply with the latest issue of SANS 149; Compliance certificate to be submitted with bid. | | |
| The following must be noted on the packaging and box: Trade name; Size and specification; Manufacturing site; CE number; Lot number; Expiry date | | |
| Packaging: Pouch of 5 balls; Unit = 1 pouch | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 568 05 | Cleansing ball swabs - X-Ray detectable, Unsterile, Diameter: 65mm Packaging: Box of 1000 balls |
| | 30 568 06 | Cleansing ball swabs - X-Ray detectable, Sterile; Diameter: 65mm 5 balls in an easy to open peel pouch Packaging: 1 unit = 50 Pouches |

**PHARMACY
CAPSULE BAGS**

RE-SEALABLE TABLETS/CAPSULES BAGS

Re-sealable dispensing Tablets/Capsules Bags

Thickness \geq 60 μ m

Overprinted conventional format to state the following information with a writable panel that allows for recording important information such as the dispenser, patient, medication, and dosage.

For State use only; Store in a cool dry place below 25°

KEEP OUT OF REACH OF CHILDREN/ BEKA KUDE EBANTWANENI/ HOU BUIE BEREIK VAN KINDERS

Name: _____; Ref No: _____; Date: _____; EXP: _____ Batch No: _____

Medicine: _____

Thatha/Take : Tablets/Capsules/ Iphilisi ; Ngelenga/ Times a day ; Every/ Emva Hour/ Ihora

Take: BEFORE/ WITH/ AFTER MEALS; PHAMBI/ NOKUDLA/ EMVA KOKUDLA

Symbols for: Sunrise/ Daytime/ Sunset/ Evening

Each facility will require additional information to be printed on the bags: **Facility Name and Address; Telephone number.**

Packaging: 100 resealable bags; Unit = 1 bale of 100 bags

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 74 010 05 | Re-sealable Tablets/Capsules Bags Size: 1 Length: \pm 100mm Width: \pm 104mm Lip: \pm 10mm Overprinted conventional format |
| | 74 010 15 | Re-sealable Tablets/Capsules Bags Size: 2 Length: \pm 120mm Width: \pm 104mm Lip: \pm 10mm Overprinted conventional format |
| | 74 010 45 | Re-sealable Tablets/Capsules Bags Size: 3 Length: \pm 160mm Width: \pm 104mm Lip: \pm 10mm Overprinted conventional format |

TABLET VIALS

| TABLET VIAL WITH HINGED CLOSURE - Clear Amber | | |
|--|--------------|---|
| Manufactured from clear amber polypropylene . | | |
| The containers must conform to the latest U.S.P. (unique selling point) specifications for tightly closed containers in respect of light transmission and vapour permeability. | | |
| Packaging: Unit = 1 box of 100 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 74 200 11 | Tablet Vial with hinged closure 10ml |
| | 74 200 21 | Tablet Vial with hinged closure 20ml |
| | 74 200 31 | Tablet Vial with hinged closure 30ml |
| | 74 200 51 | Tablet Vial with hinged closure 50ml |
| | 74 200 61 | Tablet Vial with hinged closure 75ml |
| | 74 200 71 | Tablet Vial with hinged closure 100ml |

BROWN PAPER BAG

| BROWN PAPER BAGS | | |
|---|--------------|--|
| Bags Brown paper | | |
| To comply with the latest issue of SANS 1347:2014 and proof of compliance must be submitted. | | |
| Dimensions allowed with a 15% variance. | | |
| Packaging: Bale of 100 bags; Unit = 1 bale of 100 bags | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 40 017 40 | Bags Brown paper Size: 1/2 Width: 75mm ; Gusset: 45mm ; Length: 150mm |
| | 40 017 42 | Bags Brown paper Size: 1 Width: 80mm ; Gusset: 60mm ; Length: 180mm |
| | 40 017 46 | Bags Brown paper Size: 2 Width: 95mm ; Gusset: 60mm ; Length: 110mm |
| | 40 017 48 | Bags Brown paper Size: 3 Width: 128mm ; Gusset: 74mm ; Length: 195mm |
| | 40 017 50 | Bags Brown paper Size: 4 Width: 124mm ; Gusset: 74mm ; Length: 225mm |
| | 40 017 58 | Bags Brown paper Size: 6 Width: 155mm ; Gusset: 90mm ; Length: 301mm |
| | 40 017 60 | Bags Brown paper Size: 8 Width: 165mm ; Gusset: 102mm ; Length: 301mm |
| | 40 017 62 | Bags Brown paper Size: 10 Width: 165mm ; Gusset: 102mm ; Length: 345mm |
| | 40 017 66 | Bags Brown paper Size: 12 Width: 180mm ; Gusset: 102mm ; Length: 335mm |
| | 40 017 68 | Bags Brown paper Size: 16 Width: 196mm ; Gusset: 115mm ; Length: 385mm |
| | 40 017 70 | Bags Brown paper Size: 20 Width: 210mm ; Gusset: 127mm ; Length: 390mm |
| | 40 017 72 | Bags Brown paper Size: 25 Width: 208mm ; Gusset: 127mm ; Length: 430mm |

PLASTIC BOTTLES

| PVC MEDICAL BOTTLE WITH SCREW CAP | | |
|--|-----------|--|
| Medical bottle manufactured from PVC with a screw cap made from Polypropylene . Packaging: Unit = 1 box of 100 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 74 110 22 | Plastic Rectangular Medical Bottle with screw cap Capacity: 50ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 74 110 41 | Plastic Rectangular Medical Bottle with screw cap Capacity: 100ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 74 110 85 | Plastic Rectangular Medical Bottle with screw cap Capacity: 200ml Colour: Amber Internal neck diameter 16.5mm Screw cap: 24mm |
| | 74 110 97 | Plastic Round Medical Bottle with screw cap Capacity: 500ml Colour: Amber Screw cap: 28mm Screw neck 28mm |
| | 74 110 98 | Plastic Round Medical Bottle with Spray Screw Cap Capacity: 500ml Colour: Amber |
| | 74 110 99 | Plastic Round Medical Bottle with Pump Screw Cap Capacity: 500ml Colour: Amber |
| | 74 130 30 | Plastic Round Bottle with Spray Screw Cap Capacity: 100ml Colour: Opaque White |

| POLYETHYLENE BOTTLE WITH SCREW CAP | | |
|--|-----------|---|
| Jerry Can manufactured from high density POLYETHYLENE with a screw cap made from Polypropylene . With a handle for easy carrying Packaging: Unit = Pack of 10 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 74 120 21 | Polyethylene bottle with handles and a screw cap Capacity: 2.5l Colour: Opaque White |

GLASS BOTTLES

| BOTTLE WITH SQUEEZE DROPPER | | |
|---|-----------|--|
| Medical Bottle with a dropping pipette screw-top lid Colour: Amber The squeeze-type dropper must be manufactured from polyethylene. Must be autoclavable. Packaging: Box of 100; Unit = 1 box of 100 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 74 020 10 | Glass bottle with a dropping pipette screw-top lid Capacity: 10-15ml |

| GLASS MEDICAL BOTTLE WITH SCREW CAP - AMBER | | |
|---|-----------|---|
| Medical Bottle with screw cap Colour: Amber Bottle manufactured from amber glass with a plastic/polypropylene screw top lid. Packaging: Unit = 1 box of 50 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 74 050 21 | Glass bottle with a plastic screw-top lid Capacity: 50ml |
| | 74 050 31 | Glass bottle with a plastic screw-top lid Capacity: 100ml |
| | 74 050 41 | Glass bottle with a plastic screw-top lid Capacity: 200ml |
| | 74 050 51 | Glass bottle with a plastic screw-top lid Capacity: 500ml |

MEDICINE MEASURES

SQUEEZE-TYPE BOTTLE

Squeeze-type bottle made from Linear low-density polyethylene (LLDPE) with a screw cap.

Packaging: Box of 100; Unit = 1 box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 74 020 42 | Squeeze-type dropper bottle with a screw cap Capacity: 15 ml |

MEDICINE MEASURE SPOON

Manufactured from medical plastic. **5ml**; Graduated at **2.5ml and 5ml**.

Packaging: Box of 100; Unit = 1 box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 76 716 11 | Medicine Measure Spoon 5ml Graduated at 2.5ml and 5ml |

MEDICINE TUMBLER

Made from transparent medical plastic.

Graduated from **5ml to 50ml**.

Packaging: Box of 100; Unit = 1 box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 372 01 | Measuring Tumbler Graduated: 5 - 50ml |

OINTMENT CONTAINER

OINTMENT CONTAINER

Made from high density polyethylene with a press-on lid/ screw on cap.

Packaging: Unit = 1 box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 74 260 29 | Ointment Container with snap-on Cap Capacity: 30g Colour: Clear |
| | 74 260 30 | Ointment Container with snap-on Cap Capacity: 50g Colour: White Opaque |
| | 74 260 31 | Ointment Container with snap-on Cap Capacity: 100g Colour: White Opaque |
| | 74 260 45 | Ointment Container with screw on lid Capacity: 500ml Colour: White Opaque |

SUNDRIES

| | | |
|--|--|--|
| SPUTUM CUP WITH LID | | |
| Disposable cup made from medical plastic. Capacity: approximately 150ml Packaging: Box of 100; Unit = 1 box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--------------------------------------|
| | 30 208 03 | Sputum Cups with lid. ± 150ml |

| | | |
|--|--|--|
| DENTURE CONTAINER | | |
| Disposable Round translucent plastic container with a plastic moulded lid Capacity: approximately 250ml Packaging: Box of 100; Unit = 1 box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|----------------------------------|
| | 30 092 10 | Denture container ± 250ml |

| | | |
|---|--|--|
| URINE SPECIMEN COLLECTION BOTTLE | | |
| Made from transparent and rigid medical plastic with a screw on lid. Non-graduated with a capacity of ± 90ml. Each specimen bottle must come in an individual sterile peel pack. Sterile. Single use Packaging must state: Manufacture date; Expiry date; CE Number Packaging: Box of 100; Unit = 1 box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 35 025 10 | Urine Specimen collection bottle ± 90ml |

| | | |
|---|--|--|
| TONGUE DEPRESSOR | | |
| Made from hardwood - Plastic is not acceptable. Must have rounded edges on either side. Sides and surfaces must be smooth but not slippery. Autoclavable; Single use Packaging: Box of 100; Unit = 1 box of 100 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 045 70 | Tongue Depressor Width: ± 17mm ; Length: ± 150mm Thickness: ± 1.5mm |

| | | |
|---|--|--|
| ROUND WOOD STICK APPLICATOR | | |
| Made from round smooth hardwood. Autoclavable Packaging: Box of 50; Unit = 1 box of 50 | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 010 12 | Round Wood Stick Applicator Diameter: 2mm ; Length: 150mm |

STORAGE

| PLASTIC STORAGE BIN | | |
|---|--|--|
| <p>Manufactured from lightweight, robust, heavy-duty, waterproof, rustproof, oil and chemical resistant polypropylene. Must be stackable and able to handle weight load as per size stated. Variation of 15% allowed in dimension. Front façade must allow access to content and space for a label. Must be available in different colours – red; blue yellow and in series for all the sizes. State brand name</p> <p>Packaging: Each</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 74 800 21 | Plastic storage bin Size 1 Length: 135mm Width: 105mm Height: 70mm For loads \leq 3kg |
| | 74 800 22 | Plastic storage bin Size 2 Length: 185mm Width: 105mm Height: 70mm For loads \leq 3kg |
| | 74 800 23 | Plastic storage bin Size 3 Length: 275mm Width: 140mm Height: 120mm For loads \leq 10kg |
| | 74 800 24 | Plastic storage bin Size 4 Length: 380mm Width: 205mm Height: 170mm For loads \leq 15kg |
| | 74 800 25 | Plastic storage bin Size 5 Length: 380mm Width: 420mm Height: 170mm For loads \leq 20kg |

| Corrugated Cardboard Carton | | |
|---|--|--|
| <p>Purpose: Used in the pharmacy for storage and packaging of drugs Die-cut box , waxed on both sides; Off-White Corrugated board outer, B flute inner, Manual self-locking assembly, Stackable after assembly Similar or equal to Brazier Bin</p> <p>Packaging: Each Unit = 1 carton</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 39 010 15 | Corrugated Cardboard Carton Length: 305mm Width: 102mm Height: 120mm |
| | 39 010 16 | Corrugated Cardboard Carton Length: 305mm Width: 152mm Height: 120mm |
| | 39 010 17 | Corrugated Cardboard Carton Length: 300mm Width: 120mm Height: 110mm |
| | 39 010 18 | Corrugated Cardboard Carton Length: 390mm Width: 102mm Height: 120mm |

SURGICAL SUNDRIES

MEDICAL RAZOR

Purpose: For surgical shaving
Sharp, double edged with carbon steel or medical steel blade.
For single patient use.

Packaging: Unit = 1 Box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---------------|
| | 30 465 10 | Medical Razor |

STITCH CUTTER

Manufactured from stainless steel with curved and pointed cutting edges at one end.
Sterile, single use
Packed in sterile peel-apart pack – that peels evenly without tearing paper lining.

Packaging: Unit = 1 Box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 248 38 | Stitch Cutter – disposable Length: 9cm - 12cm |

BLADES and HANDLES

| |
|--|
| <p>SURGICAL BLADE</p> <p>Made from high quality carbon steel and manufactured in accordance with BS 2982 and EN 27740</p> <p>Must be stored at room temperature and be ready to use out of the package.</p> <p>Packaging: Box of 100; individually foil wrapped in a sterile easy to open packaging</p> <p>The following must be noted on the packaging:</p> <p>Size and specification; Trade Name; CE Number; Method of sterilisation; Manufacturing site; Lot Number; Manufacture Date; Expiry date</p> |
|--|

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 045 25 | Blades Surgical Scalpel Bp No 10 - for making large straight incisions through skin and subcutaneous tissues |
| | 30 045 26 | Blades Surgical Scalpel Bp No 11 - elongated triangular blade sharpened along the hypotenuse edge |
| | 30 045 27 | Blades Surgical Scalpel Bp No 12 - for suture removal - curved and pointed cutting edge at one end |
| | 30 045 28 | Blades Surgical Scalpel Bp No 15 - with a small, curved cutting edge |
| | 30 045 29 | Blades Surgical Scalpel Bp No 20 - with a curved cutting edge and a flat, unsharpened back edge |
| | 30 045 30 | Blades Surgical Scalpel Bp No 21 - with a large, curved blade |
| | 30 045 31 | Blades Surgical Scalpel Bp No 22 - with a curved cutting edge and a flat, unsharpened back edge |
| | 30 045 32 | Blades Surgical Scalpel Bp No 23 - "leaf-shaped" blade sharpened along its leading edge. |
| | 32 045 33 | Blades Surgical Scalpel Bp No 24 - semi-circular in shape sharpened along its leading edge |

| |
|---|
| <p>SURGICAL SCALPEL HANDLE</p> <p>Purpose: Designed to hold and control surgical blades during various medical procedures. Manufactured from autoclavable stainless steel. or other durable materials.</p> <p>The handle consists of two main parts: the blade attachment mechanism and the handle grip. The blade attachment mechanism must secure the scalpel while the handle grip must be comfortable and ergonomic for precise control and manoeuvrability during surgeries.</p> <p>Packaging: Per Handle</p> |
|---|

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|---------|---|
| | 60 0039 | Surgical Scalpel handle No.3 for blade sizes 10 - 17 |
| | 60 0040 | Surgical Scalpel handle No.4 for blade sizes >18 |

SKIN GRAFT

| SKIN GRAFT BLADES - HUMBY | | |
|---|-----------|--|
| <p>Purpose: For the surgical removal of an area of healthy skin from one zone of the patient's body, to be transplanted to another area Must be compatible with Humby, Watson Braithwait and Corbett knives - with 3 Holes for securement. Made from Stainless Steel that does not bend or warp during use. The sharp blade edge must allow smooth shaving of healthy skin tissue. Sterile and individually wrapped.</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 045 15 | Skin Graft Blades: Bodenham Humby 157mm x 18mm x 0.23mm |

| MESH DERMACARRIER | | |
|--|-----------|--|
| <p>Purpose: Carrier board for the meshing and expansion of harvested skin Must consist of a highly flexible medical grade polypropylene. Plateau grooves must help ensure correct alignment in the mesher. For single use; Individually packed and sterilised with ETO Gas</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 376 03 | Mesh Derma-carrier Size 1.5:1 Length: 200mm Width: 80mm |
| | 30 376 02 | Mesh Derma-carrier Size 3:1 Length: 200mm Width: 80mm |

ENDOSCOPY

| LAPAROSCOPY PORTS | | |
|---|----------|---|
| <p>Purpose: To gain access into the abdomen for visualisation and manipulation of intra-corporeal structures Must consist of a stability sleeve and insertion trocar. The sleeve must have a double valve system that prevents leakage of CO2 during the procedure. The sleeve must be ribbed externally to allow it to be held in place by the abdominal wall. The trocar must insert into the sleeve seamlessly and clip into place. The insertion trocar must have a bladeless, transparent, optical tip to allow for visualisation of abdominal wall layers and visceral contents during insertion. There must be a mechanism to hold an appropriate size lens in place during insertion.</p> <p>The following must be noted on the packaging: Size and specification; Trade Name; CE Number; Method of sterilisation; Manufacturing site; Lot Number; Manufacture Date; Expiry date Packaging: Each</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | LapPort1 | Laparoscopic Port Set - 4 Ports: Optical bladeless tip trocars and stability sleeves . 2 Ports: Length: 100mm Diameter: 5mm ; 2 Ports: Length: 100mm Diameter: 12mm |
| | LapPort2 | Laparoscopic Port - Optical bladeless tip trocar and stability sleeve . Length: 100mm Diameter: 5mm |
| | LapPort3 | Laparoscopic Port - Optical bladeless tip trocar and stability sleeve . Length: 100mm Diameter: 10mm |

STAPLERS

LINEAR CUTTER STAPLER

Purpose: Gastrointestinal anastomosis stapler, used for resection and transection of organs or tissues

Consists of disposable Linear Cutter Stapler that is loaded with two double staggered rows of titanium staples that simultaneously cut and divide tissue between the two double rows linear cutter.

Must be easily to operate by squeezing of the handles followed by shifting of side knobs forth and back.

Must have a precise mechanism to facilitate parallel jaws closing followed by proper staples formation.

Must be reloadable at least 7 times to allow for 8 firings per instrument.

Each new cartridge should have a new blade for precise incision.

The following must be noted on the packaging:

Size and specification; Trade Name; CE Number; Method of sterilisation; Manufacturing site; Lot Number; Manufacture Date; Expiry date

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|----------|---|
| | LCStap | Linear Cutter Stapler 80mm Width: 3.8mm |
| | LCReload | Linear Cutter Stapler Reloads 80mm Width: 3.8mm |

SKIN STAPLES GUN

Pre-Packed with **35 staples** manufactured from surgical steel Width: **5.6mm**; Height: **3.8mm**

Must offer a clear view of the operating site and it must be easy to gauge the number of remaining staples.

Staples must **approximate completely**. There must be audible feedback when the staple is deployed.

Sterile and individually packed in a peel pouch that is easy to open.

For single use only

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Box of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 473 26 | Skin Staples Gun – Single use Pre-Packed with 35 staples Width: 5.6mm ; Height: 3.8mm |

SKIN STAPLES REMOVER - DISPOSABLE

Sterile and individually packed in a peel pouch that is easy to open.

For single use only

The following must be noted on the packaging:

Trade name; Manufacturing site; Expiry date; Method of sterilisation

Packaging: Box of 10; Unit = 1 box of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 473 24 | Skin Staples Remover - DISPOSABLE |

SKIN STAPLES REMOVER - REUSABLE

Manufactured from Surgical steel - Autoclavable.

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 473 25 | Skin Staples Remover - REUSABLE |

DIATHERMY

DIATHERMY PLATES

Purpose: For use in electrosurgery applications.

Consists of a highly conductive adhesive grounding pad that forms part of the circuit allowing high density radio frequency current to pass through tissue and to cut or coagulate tissue. Must be suitable for all types of surgery.

Must have a large conductive surface and a 3m 3-pin cable.

There must be a pull tab to help removal of the from the patient.

Single use, individually packaged, hypoallergenic, with an adhesive border to help to prevent fluid intrusion.

The following must be noted on the packaging:

Size and specification; Trade Name; CE Number; Manufacture Date; Expiry date

Packaging: 1 unit = Pack of **50 plates**

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|-------------------------------------|
| | 30 249 98 | Diathermy Plates Child D135P |
| | 30 249 99 | Diathermy Plates Adult |

DIATHERMY HAND CONTROL

Purpose: For electrosurgical cutting and haemostasis during medical procedures

Consists of a slim, lightweight, and ergonomic electrosurgical Pencil that forms part of the diathermy circuit.

It must be connected to the energy driver via a **3m, 3-pin cable**.

The pencil must have a sealed design that protects the inside circuit board.

It must have push-button switches that allow for good tactile feedback - preventing accidental activation.

Surgeons should be able to hear and feel when the switch engages.

Must allow for accurate cutting to make the tissue resection cleaner and the incision smooth.

Single-use and packaged in sterile easy to open packaging.

The following must be noted on the packaging:

Size and specification; Trade Name; CE Number; Manufacture Date; Expiry date

Packaging: 1 unit = Pack of **50**

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|------------------------|
| | 30 250 14 | Diathermy Hand Control |

DIATHERMY SCRATCH PADS

Purpose: For cleaning diathermy knives

Consists of an abrasive pad with an adhesive backing

Sterile and individually packed in a peel pouch that is easy to open.

Packaging: Box of 100; Unit = 1 box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 440 30 | Diathermy Scratch Pads Width: 8cm Length: 10cm |

FEEDING

NASOJEJUNAL FEEDING TUBES

Purpose: Naso-jejunal feeding

The feeding tube must be soft, kink resistant, **radiopaque** and graduated with centimetre markings.

Must have an open distal end with 2 lateral eyes and atraumatic tip.

Must have an **introducer/ stylet** and must be **weighted** end for easy placement.

The proximal end must fit a luer syringe and have a plug to close the end.

To comply with **ISO/DIS 20695:2017** or equivalent.

Manufactured from latex free non- toxic and non- irritant **silicone/polyurethane** and medical grade plastic.

Must be sterile and individually packed in a peel pouch that is easy to open; For single use only.

To comply with ISO/DIS 20695:2017 or equivalent.

The following must be noted on the packaging:

Trade name ; Size and specification; Method of sterilisation; Manufacturing site ; CE number; Lot number; Expiry date

Packaging: Unit = 1 box of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 33 006 69 | Naso-jejunal feeding tube - 6Fr Length: 550mm +/- 10% |
| | 33 006 70 | Naso-jejunal feeding tube - 8Fr Length: 550mm +/- 10% |
| | 33 006 71 | Naso-jejunal feeding tube – 10Fr Length: 1100mm +/- 10% |
| | 33 006 72 | Naso-jejunal feeding tube – 12Fr Length: 1100mm +/- 10% |

NASOGASTRIC TUBE - RYLES TYPE

Purpose: For aspiration of stomach contents or introduction of feeds into the stomach

The Ryles Nasogastric Tube must be soft yet kink resistant, must remain firm at body temperature, and must not occlude when suction is applied at the proximal end.

Must have a radiopaque line and be graduated at 45/50, 55/60, 65/70, and/or 75cm from the distal tip for ease of placement.

Must have a closed, smooth, and moulded distal tip with 4 lateral eyes.

The distal tip must be weighted and coned with corrosion resistant, radio-opaque stainless-steel balls sealed into the tube, or it must have a smooth moulded stainless-steel weight.

Must have a proximal end that is colour coded and funnel shaped with a luer adapter and a plug to close the proximal end.

The proximal funnel must be firm and must not collapse without the luer adaptor and plug in place.

Must be manufactured from latex free non- toxic and non- irritant transparent/frosted PVC and medical grade plastic.

Must be sterile and individually packed in a peel pouch that is easy to open, For single use only.

To comply with ISO/DIS 20695:2017 or equivalent.

The following must be noted on the packaging:

Trade name ; Size and specification; Method of sterilisation; Manufacturing site ; CE number; Lot number; Expiry date

Packaging: Unit = 1 box of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 33 006 79 | Nasogastric tube - Ryles type – 8Fr Length: 100 cm +/- 10% Colour: Blue |
| | 33 006 80 | Nasogastric tube - Ryles type – 10Fr Length: 100 cm +/- 10% Colour: Black |
| | 33 006 81 | Nasogastric tube - Ryles type – 12Fr Length: 100 cm +/- 10% Colour: White |
| | 33 006 82 | Nasogastric tube - Ryles type – 14Fr Length: 100 cm +/- 10% Colour: Green |
| | 33 006 83 | Nasogastric tube - Ryles type – 16Fr Length: 100 cm +/- 10% Colour: Orange |
| | 33 006 84 | Nasogastric tube - Ryles type – 18Fr Length: 100 cm +/- 10% Colour: Red |
| | 33 006 85 | Nasogastric tube - Ryles type – 20Fr Length: 100 cm +/- 10% Colour: Yellow |

GASTRIC

FLATUS TUBE

Purpose: To decompress the colon

Made from medical grade PVC.

Distal end: Open end with ONE lateral eye **Proximal end:** Funnel connector.

Must be sterile and individually packed in a peel pouch that is easy to open, For single use only.

The following must be noted on the packaging:

Trade name ; Size and specification; Method of sterilisation; Manufacturing site ; CE number; Lot number; Expiry date

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 33 135 17 | Flatus Tube 28Fr Minimum length: 500mm |

STOMACH WASH OUT TUBE

Purpose: To decompress the stomach and allow for washout of large matter

Made from medical grade, non-toxic PVC.

Distal end: Rounded tip with four eyes.

Must be sterile and individually packed in a peel pouch that is easy to open, For single use only.

The following must be noted on the packaging:

Trade name ; Size and specification; Method of sterilisation; Manufacturing site ; CE number; Lot number; Expiry date

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 33 150 13 | Stomach Wash Out Tube 28Fr Minimum Length: 800mm |
| | 33 150 17 | Stomach Wash Out Tube 30Fr Minimum Length: 800mm |
| | 33 150 19 | Stomach Wash Out Tube 32Fr Minimum Length: 800mm |

SENGSTAKEN TUBE

Purpose: Tamponade of gastro-oesophageal bleeding that is unresponsive to medical and endoscopic treatment

Must have separate gastric and oesophageal balloons, with individual inflation ports for each balloon.

The oesophageal balloon port must inflate a small balloon in the oesophagus, and the gastric balloon port must inflate a balloon in the stomach. The gastric aspiration port must remove fluid and air out of the stomach.

Latex free, DEHP free.

Must be sterile and individually packed in a peel pouch that is easy to open, For single use only.

The following must be noted on the packaging:

Trade name ; Size and specification; Method of sterilisation; Manufacturing site ; CE number; Lot number; Expiry date

Packaging: Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | Sengs | Sengstaken Tube 21Fr Length: 100cm |

WOUND DRAINS

CORRUGATED DRAINS

Purpose: Corrugated sheets that can be cut to the required size for the passive drainage of blood, serum, pus, or urine.

Manufactured from latex and DEHP free silicone elastomer.
 Must not contain substances that can be extracted by water or tissue fluids.
 Must be water- and tissue-repellent, reducing adhesion to blood or tissue matter.
 Sterile - for single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Each; Unit = 1 drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 226 10 | Corrugated Drains - 350mm x 45mm (±5%) |
| | 30 226 11 | Corrugated Drains - 350mm x 110mm (±5%) |

CORRUGATED DRAIN WITH LUMENS

Purpose: Corrugated sheets with internal lumens that can be cut to the required size for the passive drainage of blood, serum, pus, or urine.

Manufactured from latex- free and DEHP free silicone elastomer with lumens of **3mm**; 12 Joined catheter lumens per sheet.
 Must not contain substances that can be extracted by water or tissue fluids.
 Must be water- and tissue-repellent, reducing adhesion to blood or tissue matter.
 Must be Single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Each; Unit = 1 drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 226 12 | Corrugated Drain Sheet with drainage lumens Size: 40cm x 110mm (±5%) |

PENCIL DRAINS

Purpose: Corrugated tube for the passive drainage of blood, serum, pus, or urine from smaller areas

Manufactured from latex and DEHP free silicone elastomer.
 Must not contain substances that can be extracted by water or tissue fluids.
 Must be water- and tissue-repellent, reducing adhesion to blood or tissue matter.
 Sterile - for single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Each; Unit = 1 drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|-------------------------------------|
| | PD1 | Pencil Drains - 6mm x 300mm |
| | PD2 | Pencil Drains - 10mm x 300mm |
| | PD3 | Pencil Drains - 20mm x 300mm |

HUBLESS WOUND DRAIN WITH TROCAR

Purpose: A radiopaque drain used to remove excess fluid from a cavity.

Manufactured from a soft flexible white, radiopaque silicone with a solid centre and four channels running down the sides. Must have a blue radiopaque stripe running along the length of the drain, round silicone extension tube and adapter. Must have a trocar attached.

Must be latex and DEHP free; Sterile - for single use - packed in sterile peel-apart pack.

Single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Per Drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | WDR1 | Wound Drain - Hubless 10 FR Round, 1/8" Trocar |
| | WDR2 | Wound Drain - Hubless 15 FR Round, 3/16" Trocar |
| | WDR3 | Wound Drain - Hubless 19 FR Round, 1/4" Trocar |
| | WDR4 | Wound Drain - Hubless 19 FR Round, 1/4" Bendable Trocar |

HUBLESS WOUND DRAIN WITHOUT TROCAR

Purpose: A radiopaque drain used to remove excess fluid from a cavity.

Manufactured from a soft flexible white, radiopaque silicone with a solid centre and four channels running down the sides. Must have a blue radiopaque stripe running along the length of the drain, round silicone extension tube and adapter.

Must be latex and DEHP free; Sterile - for single use - packed in sterile peel-apart pack.

Single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Per Drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | WDR9 | Wound Drain - Hubless, 24Fr (8.0mm) Round |

FLAT FLUTED WOUND DRAIN

Purpose: A radiopaque drain used to remove excess fluid from a cavity.

Manufactured from a soft flexible white, radiopaque silicone with a solid centre and four channels running down the sides. Must have a blue radiopaque stripe running along the length of the drain, round silicone extension tube and adapter.

Must be latex and DEHP free; Sterile - for single use - packed in sterile peel-apart pack.

Single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Per Drain

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | WDF1 | Wound Drain – 10mm Flat, 3/4 Fluted |
| | WDF2 | Wound Drain – 10mm Flat, 3/4 Fluted with 3/16" with a Trocar |
| | WDF3 | Wound Drain – 10mm Flat, Full Fluted |
| | WDF4 | Wound Drain – 10mm Flat, Full Fluted with 3/16" with a Trocar |

RESERVOIRS

CLOSED WOUND DRAIN RESERVOIR

Purpose: Reservoir with a clear window to view fluid being drained from a surgical wound site

Consists of a closed reservoir with an anti-reflux valve and Y -connector and is compatible with all wound drains.
Sterile - for single use - packed in sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Unit = pack of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | CWDR1 | Closed wound drainage reservoir 150ml |
| | CWDR2 | Closed wound drainage reservoir 300ml |
| | CWDR3 | Closed wound drainage reservoir 450 - 500ml |

BULB SUCTION RESERVOIR

Purpose: Clear reservoir with a clear view of fluid being suctioned from a surgical wound site

100ml bulb reservoir closed reservoir with an **anti-reflux valve** Must be compatible with **“Wound drains”**.
Must be sterile for single use and packed in a sterile peel-apart pack.

The following must be noted on the packaging:

Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date

Packaging: Unit = pack of 10

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--------------------------------------|
| | CWDR4 | Bulb Suction Reservoir 100 ml |

OPHTHALMOLOGY

| EYE SHIELD | | |
|--|----------|--|
| Purpose: To protect operated or threatened eyes Manufactured from a clear polycarbonate material. The shield must have nine holes in its centre, have smooth surfaces and a rounded, lipped edge. There must be no string or string-like attachments. The following must be noted on the packaging: Trade name; Size and specification; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | ESheildA | Eye Shield - Adult 6.2cm x 7.5cm |
| | ESheildP | Eye Shield - Paediatric 5.1cm x 4.1cm |

| EYE CONFORMER | | |
|--|--------|--|
| Purpose: To protect operated or threatened eyes Manufactured from a clear acrylic material. The conformer must have holes in its centre, have smooth surfaces and edges. Must be curved to fit the natural curve of the eye socket. Single use, sterile and individually packed in a peel pouch that is easy to open. The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10 | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | EConf1 | Eye Conformer Size: 19mm x 22mm |
| | EConf2 | Eye Conformer Size: 21mm x 24mm |
| | EConf3 | Eye Conformer Size: 22mm x 25mm |

| ORBITAL BALL IMPLANT | | |
|---|-------|---|
| Purpose: To maintain the contour of the eyeball, maintain eye movements, provide a base for prosthesis, and maintain shape of eye socket and to prevent bony deformity. Manufactured from silicone. Must have a smooth non irritating surface. Single use, sterile and individually packed in a peel pouch that is easy to open. The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: each | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | Orb1 | Orbital Ball Implant Size: 10-14mm |
| | Orb2 | Orbital Ball Implant Size: 16mm |
| | Orb3 | Orbital Ball Implant Size: 18mm |
| | Orb4 | Orbital Ball Implant Size: 20mm |

OPHTHALMIC KNIVES

| OPHTHALMIC CRESCENT KNIFE | | |
|--|-----------|----------------------------------|
| <p>Purpose: For Ophthalmic Surgery</p> <p>Consists of a micro-crescent blade with the plane of the blade angled upward at 45 degrees to the handle. The micro-crescent blade must have cutting surfaces at the knife tip as well as on both sides. The cutting surfaces at the tip of the micro-crescent blade must be rounded, semi-rounded or squared off. The maximum width of the micro-crescent blade must be 2.2mm or less. Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | CresKnife | Ophthalmic Crescent Knife |

| OPHTHALMIC BLEB KNIFE | | |
|---|-----------|------------------------------|
| <p>Purpose: For Ophthalmic Glaucoma Surgery</p> <p>Consists of a micro- blade with a straight plane and a maximum width of 1mm. The micro-crescent blade must have cutting surfaces at the knife tip as well as on both sides. The cutting surfaces at the tip of the blade must be rounded, semi-rounded or squared off. The knife must have an anti-glare surface. Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | BlebKnife | Ophthalmic Bleb Knife |

| OPHTHALMIC SLIT KNIFE | | |
|--|-----------|------------------------------|
| <p>Purpose: For Ophthalmic Surgery</p> <p>Consists of a micro- blade with the plane of the blade angled upward at 45 degrees to the handle. The blade must have cutting surfaces at the knife tip as well as on both sides. The cutting surfaces must be anti-glare. The maximum width of the blade must be 2.75mm. Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | SlitKnife | Ophthalmic Slit Knife |

| OPHTHALMIC STAB KNIFE | | |
|--|-----------|------------------------------|
| <p>Purpose: For Ophthalmic Surgery</p> <p>Consists of a micro- blade with the plane of the blade angled at 15 degrees to the handle. The blade must have cutting surfaces at the knife tip as well as on both sides. The maximum width of the blade must be 1mm or less. Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | StabKnife | Ophthalmic Stab Knife |

| OPHTHALMIC SPEAR INCLUDING STICK | | |
|--|-------|----------------------------------|
| <p>Purpose: For quick absorption of blood & fluid from the surgical area.</p> <p>Consists of an arrow shaped spear made of a plastic stick and a triangular compressed absorbent sponge. The sterile spear must be made from highly absorbent natural cellulose material for wicking and fluid containment. The spear must have a flat and rigid tip to allow precision absorption during delicate surgical procedures and must have a fast-wicking action.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Unit = pack of 10</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | Spear | Ophthalmic Spear Including Stick |

| TRYPAN BLUE OPHTHALMIC SOLUTION | | |
|---|--------|---------------------------------|
| <p>Purpose: To stain anterior capsule during cataract surgery</p> <p>Consists of a 0.06% solution containing 0.6 mg trypan blue; 1.9 mg sodium mono-hydrogen orthophosphate (Na₂HPO₄·2H₂O); 0.3 mg sodium di-hydrogen orthophosphate (NaH₂PO₄·2H₂O); 8.2 mg sodium chloride (NaCl); and water for injection. The pH is 7.3 - 7.6. The osmolality is 257-314 mOsm/kg.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Each</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | Trypan | Trypan Blue Ophthalmic Solution |

| MINI MONOKA STENT | | |
|---|--|--|
| <p>Purpose: To treat canalicular lacerations involving the external 2/3 of one canaliculus.</p> <p>Silicone tubing Length: 40mm ID: 0.3mm OD: 0.64mm Must have a proximal punctum fixation device collar with a width of 2mm. Must be pyrogen-free and toxin-free. Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Each</p> | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|----------|--|
| | MonStent | Mini Monoka Stent Length: 40mm ID: 0.3mm OD: 0.64mm |

| GLAUCOMA DRAINAGE DEVICE - AHMED VALVE |
|---|
| <p>Purpose: Emergency treatment of glaucoma and the immediate reduction of intraocular pressure</p> <p>Consists of a valved plate body; a valve with a casing and a draining tube Plate/Valve: Thickness: 2.1mm Width: 13mm Length: 16mm; Surface Area: 184mm² Drainage Tube: Length: 25.4mm ID: 0.305mm OD: 0.635mm Must have a tapered profile for easy insertion as a single-stage procedure. Must not obstruct, cause excessive drainage, or result in chamber collapse. Made of medical grade silicone and polypropylene Single use, sterile and individually packed in a peel pouch that is easy to open.</p> <p>The following must be noted on the packaging: Trade name; Size and specification; Method of sterilisation; Manufacturing site; CE number; Lot number; Expiry date Packaging: Each</p> |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|--------|--|
| | AValve | Glaucoma Drainage Device - Ahmed Valve |

CSSD

TAPE

| AUTOCLAVE SEALING TAPE | | |
|--|--------------|---|
| Purpose: Adhesive Tape to indicate whether a specific temperature has been reached in an autoclave Must change colour after exposure to temperatures commonly used in sterilisation processes, typically 121°C in a steam autoclave. To comply with the latest issue of CKS 224 Roll | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 577 49 | Autoclave Indicator Tape Size: 12mm x 50m |
| | 30 577 50 | Autoclave Indicator Tape Size: 24mm x 50m |

| AUTOCLAVE ADHESIVE TAPE | | |
|--|--------------|--|
| Purpose: Adhesive Tape to seal Autoclave steam bags Must remain adherent during the autoclaving process. Roll | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 30 578 05 | Autoclave Adhesive tape Size: 24mm x 50m |

| Instrument Colour Code Tape | | |
|--|--------------|---|
| Purpose: Unique vinyl tape for marking all types of surgical instruments. Instrument Colour Code Tape. Must be autoclavable and withstands repeated sterilisation cycles. Must be extremely durable and can be used in all types of autoclaves | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | CCTWhite | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: White |
| | CCTBlue | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Blue |
| | CCTGreen | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Green |
| | CCTRed | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Red |
| | CCTYellow | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Yellow |
| | CCTPurple | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Purple |
| | CCTOrange | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Orange |
| | CCTBrown | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Brown |
| | CCTBlack | Instrument Colour Code Tape Size: ± 6.35mm x 10m Colour: Black |

INDICATORS

STEAM STERILISATION INDICATOR STRIPS

Purpose: Paper strip printed with a chemical indicator ink that turns from off-white to dark brown/black when exposed to the steam sterilisation process

Internal chemical indicator paper strip pack must indicate whether steam has penetrated to the point of placement of the strip, usually the centre of a pack.

Each strip is perforated at the centre if a shorter internal chemical indicator is preferred.

Pack of 240 strips

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 338 35 | Steam Sterilisation Indicator Strips Size: 1.5cm x 20cm |

BOWIE-DICK INDICATOR

Purpose: Vacuum test pack to test the efficiency of **porous load autoclaves**

The pack must be made from virgin pulp paper holding in its centre the chemical indicator sheet.

The entire pack must be strengthened by a paper wrap binding, which guarantees the constant compactness of the test pack.

Must have an external chemical indicator.

Must fully comply with the latest edition of ISO 11140-4; ISO Compliance Certificate to be submitted with bid.

Each

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 338 32 | Bowie Dick Indicator for Autoclaves |

ETHYLENE OXIDE (EO) CHEMICAL INDICATOR STRIP

Purpose: To indicate whether EO sterilisation conditions have been met at the point of placement of the strip

Consists of a paper strip printed with a chemical indicator ink that turns from red to green when exposed to the critical parameters of the EO sterilisation process (i.e., gas concentration, RH, time, and temperature).

Must indicate whether EO sterilisation conditions have been met at the point of placement of the strip, which should be placed in the area determined to be the least accessible to EO penetration.

Each strip is perforated at the centre if a shorter internal chemical indicator is preferred.

Pack of 240 strips

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 338 30 | Ethylene Oxide (EO) Chemical Indicator Strip (Attest strip) Size: 1.5cm x 20cm |

BIOLOGICAL INDICATOR TEST TUBE

Purpose: A self-contained rapid readout system for steam sterilisation processes

Consists of a polypropylene tube and cap containing a filter and purple culture medium

The culture medium must indicate if the sterilisation process has been successful - Successful sterilisation must result in the medium remaining purple.

Failure will result in the medium changing from purple to yellow to green - showing the presence of live spores.

Readout after incubation at 60°C; Must have a 2-year shelf life

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|--|
| | BioTT | Biological Indicator Test Tube (BT220 equivalent) |

BIOLOGICAL INDICATOR THERMAL PRINTING PAPER

Purpose: To record the results of biological indicator subjected to sterilisation

White thermal paper roll Diameter: **50mm**; Width **58mm**; Must have a 2-year shelf life

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-------|---|
| | BioPP | Biological Indicator Thermal Printing Paper |

STERILISATION BAGS

PAPER BAG FOR STEAM STERILISATION

Purpose: Must have a process indicator pattern/band that changes colour to indicate when the properly steam autoclaved

Manufactured from medical graded paper with a process indicator that changes colour when exposed to saturated steam .

PRINTING ON BAG:

FRONT: Do not use if the wrapper or seal is damaged, Indicator pattern/band and final colour of indicator ink.

BACK: Manufacturers name or trademark, Size of bag, Designation, Batch number

Paper to comply with ISO 11607; SANS 50868; DIN 58953-6; CKS 595 or equivalent; ISO 9001

PACKAGING: bale/box of 1 000 bags,

Packed in a suitable rigid cardboard box/bale

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 018 10 | Paper Bag for Steam Sterilisation Size: 1 Width: 90mm Length: 100mm Gusset: 50mm |
| | 30 018 11 | Paper Bag for Steam Sterilisation Size: 2 Width: 90mm Length: 230mm Gusset: 50mm |
| | 30 018 12 | Paper Bag for Steam Sterilisation Size: 3 Width: 115mm Length: 190mm Gusset: 35mm |
| | 30 018 13 | Paper Bag for Steam Sterilisation Size: 4 Width: 125mm Length: 260mm Gusset: 50mm |
| | 30 018 15 | Paper Bag for Steam Sterilisation Size: 6 Width: 190mm Length: 320mm Gusset: 60mm |
| | 30 018 16 | Paper Bag for Steam Sterilisation Size: 7 Width: 180mm Length: 380mm Gusset: 75mm |
| | 30 018 17 | Paper Bag for Steam Sterilisation Size: 8 Width: 65mm Length: 480mm Gusset: 30mm |
| | 30 018 18 | Paper Bag for Steam Sterilisation Size: 9 Width: 120mm Length: 135mm Gusset: 30mm |
| | 30 018 19 | Paper Bag for Steam Sterilisation Size: 11 Width: 290mm Length: 445mm Gusset: 125mm |
| | 30 018 20 | Paper Bag for Steam Sterilisation Size: 12 Width: 330mm Length: 560mm Gusset: 150mm |
| | 30 018 31 | Paper Bag for Steam Sterilisation Size: 1020 Width: 65mm Length: 200mm Gusset: 35mm |

STERILISATION PAPER

Crêped Sterilisation Paper

Purpose: Used as a wrapping material for sterilisation of dressings, instruments, etc

White/Blue/Green paper with a mass not less than **57g/m²**,

Must be a moisture repellent and resistant to air and fluid borne barrier.

Must comply with the latest issue of SANS 50868-2, DIN 58953-6, SANS 11607-1 and/or EN868-2 and DIN 59953-6

Box of 100

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 445 32 | Crêped Sterilisation Paper White / Blue / Green Size: 300mm x 300mm |
| | 30 445 31 | Crêped Sterilisation Paper White / Blue / Green Size: 600mm x 600mm |
| | 30 445 33 | Crêped Sterilisation Paper White / Blue / Green Size: 600mm X 900mm |
| | 30 445 36 | Crêped Sterilisation Paper White / Blue / Green Size: 900mm x 900mm |
| | 30 445 35 | Crêped Sterilisation Paper White / Blue / Green Size: 900mm x 1300mm |
| | 30 445 34 | Crêped Sterilisation Paper White / Blue / Green Size: 900mm x 1500mm |
| | 30 445 29 | Crêped Sterilisation Paper White / Blue / Green Size: 1000mm x 1400mm |
| | 30 445 27 | Crêped Sterilisation Paper White / Blue / Green Size: 1500mm x 1500mm |

STERILISATION WRAP

| STERILISATION WRAP - SINGLE REGULAR (SPUNBOUND) | | |
|---|-----------|--|
| Sterilisation Wrap - Single Regular Medical, non-woven Spunbond-Meltblown-Spunbond (SMS) To comply with the latest issue of SANS 50868-2, DIN 58953-6, SANS 11607-1 and/or EN868-2 Pack of 240 Sheets | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | 39 127 49 | Sterilisation Wrap (Spunbound) - Regular Size: 600mm x 600mm |
| | 39 127 50 | Sterilisation Wrap (Spunbound) - Regular Size: 900mm x 900mm |
| | 39 127 51 | Sterilisation Wrap (Spunbound) - Regular Size: 1200mm x 1200mm |
| | 39 127 52 | Sterilisation Wrap (Spunbound) - Regular Size: 1400mm x 1400mm |

STERILISATION REELS

| FLAT STERILISATION REELS WITHOUT GUSSET | | |
|--|--|--|
| Purpose: Suitable for packaging medical instruments, which must be sterilised in steam or ethylene oxide. | | |
| Heat sealable reels of paper and plastic film with process indicators for saturated steam and ethylene oxide. Seams to peel open without linting and/or tearing film. sterilisation process indicators for steam and ethylene oxide must be printed at equal distances throughout the roll. Product to comply with the latest issue of ISO 11607-1 and SANS 50868 Part 3 and 5. Certificates of declarations of conformance from manufacturer to be submitted: 93/42/EEC Medical Device Directive, Test report from independent agency and ISO9001 certificate of registration of quality management system. Roll | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 445 38 | Flat Sterilising Reels without gusset Size: 50mm x 200m |
| | 30 445 39 | Flat Sterilising Reels without gusset Size: 75mm x 200m |
| | 30 445 40 | Flat Sterilising Reels without gusset Size: 100mm x 200m |
| | 30 445 41 | Flat Sterilising Reels without gusset Size: 150mm x 200m |
| | 30 445 42 | Flat Sterilising Reels without gusset Size: 200mm x 200m |
| | 30 445 43 | Flat Sterilising Reels without gusset Size: 250mm x 200m |

| FLAT STERILISATION REELS WITH GUSSET | | |
|---|--|--|
| Purpose: Suitable for packaging medical instruments, which must be sterilised in steam or ethylene oxide. | | |
| Heat sealable reels of paper and plastic film with process indicators for saturated steam and ethylene oxide. Seams to peel open without linting and/or tearing film. sterilisation process indicators for steam and ethylene oxide must be printed at equal distances throughout the roll. Suitable for packaging medical instruments, which must be sterilised in steam or ethylene oxide. Product to comply with the latest issue of SANS ISO 11607-1 and SANS 50868 Part 3 and 5. Certificates of declarations of conformance from manufacturer to be submitted: 93/42/EEC Medical Device Directive, Test report from independent agency and ISO9001 certificate of registration of quality management system. Roll | | |

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|---|
| | 30 445 45 | Flat sterilising reels with gusset Size: 100mm x 150mm |
| | 30 445 46 | Flat sterilising reels with gusset Size: 100mm x 200mm |
| | 30 445 47 | Flat sterilising reels with gusset Size: 100mm x 250mm |

DETERGENTS & CLEANERS

DETERGENT LIQUID CLEANER FOR THEATRE INSTRUMENTS

Purpose: An enzymatic cleaner that removes blood, protein, mucous, vomit and faecal matter.

Must have Protease and Amylase enzymes with buffers and surfactants. (Endozime equivalent)

Must be able to access spaces and break down organic material.

Must leave instruments free of any unpleasant odours and baked-on blood.

Must work within 2 minutes

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 40 030 57 | Detergent Liquid Cleaner for Theatre instruments (Endozime equivalent) |

ENZYMATIC PRESOAK AND CLEANER FOR THEATRE INSTRUMENTS

Purpose: Used to remove a broad range of contaminations from surgical instruments and medical devices in the manual and automated cleaning process.

Enzymatic Presoak and Cleaner - must be low foaming with a powerful enzyme that shortens the length of cleaning processes.

A **Neutral Detergent** must have surfactants to penetrate all spaces and must work quickly to remove the most proteins and lipids, quickly and effectively.

Must have chelating and sequestering agents to protect instrument surfaces.

An **Alkaline Detergent** must be compatible with stainless steel, plastics, and soft metals and must also be safe for use on aluminium container systems.

Must protect the instruments from corrosion and remove hard water scale deposits, rust and stains on both instruments and the washing chamber.

Must be biodegradable, environmentally friendly.

Pack: 3.8 Litres

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|-----------|--|
| | 30 073 10 | Enzymatic Presoak & Cleaner - 3.8litres (Prolystica Equivalent) |
| | 30 073 11 | Neutral Detergent - 3.8litres |
| | 30 073 12 | Alkaline Detergent- 3.8litres |

DISINFECTANT STERILISING SOLUTION

Purpose: A high-level disinfectant for the soaking of heat-sensitive semi-critical medical devices such as endoscopic, respiratory therapy, and anaesthesia equipment.

Must achieve high-level disinfection within **12 minutes at 20 °C (68 °F)**. Cidex equivalent

Must be fast acting with long lasting efficacy: reusable for up to 14 days when monitored with CIDEX® OPA Solution Test Strips

Must not require activation or mixing.

Must be Glutaraldehyde-free

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|----------|--|
| | 60 00 09 | Disinfectant sterilising Solution (Cidex equivalent) |

SURGICAL INSTRUMENT RUST AND SCALE REMOVER

Purpose: Remove rust, spots, and stains from surgical instruments.

Must safely remove rust, stains, spotting, hard water scale and mineral deposits frequently encountered from sterilisation.

Must loosen stiff joints and locks.

Must work within **10 minutes**

| ICN NO: | ITEM: | DESCRIPTION: |
|---------|---------|--|
| | SurgiSt | Surgical Instrument rust and scale remover (Surgistain equivalent) |

| INSTRUMENT BRUSH | | |
|--|---------|---|
| <p>Purpose: To remove debris and contaminants, ensuring proper instrument function and reducing the risk of cross-contamination.</p> <p>Made from high grade Nylon for instruments or Stainless Steel for the effective cleaning of instrument boxes and serrations. Must have at least three rows of bristles that effectively remove debris, contaminants, and biofilm, ensuring optimal performance and patient safety in healthcare settings.</p> <p>Must be autoclavable</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | NBrush | Nylon Instrument Brush |
| | SSBrush | Stainless Steel Instrument Brush |

CSSD PROTECTION

| HEAT RESISTANT APRON | | |
|--|---------|-----------------------------|
| <p>Purpose: Protective Heat-resistant Apron</p> <p>Must consist of flexible; wipeable fabric that is resistant to radiated heat up to 1000°C</p> <p>Must secure and be adjustable around the neck and waist and be comfortable and well fitting.</p> <p>Must be at least 5mm thick</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | HRApron | Heat-resistant Apron |

| HEAT RESISTANT AUTOCLAVE GLOVES | | |
|---|----------|--|
| <p>Purpose: Protective Heat-resistant Autoclave Gloves for the handling of hot objects from autoclaves or ovens.</p> <p>Must be heat-resistant up to 232 °C</p> <p>The gauntlets must protect the wrists and forearm and must have a comfortable fit.</p> <p>Must be machine washable.</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | HRGloves | Heat-resistant Autoclave Gloves |

| NITRILE DECONTAMINATION GLOVES | | |
|---|-----------|--|
| <p>Purpose: Decontamination gloves that offer extra protection in areas where there is a risk of exposure to sharps and chemicals</p> <p>Must be ambidextrous; powder-free, non-sterile Nitrile gloves.</p> <p>Must have full finger texture and have 40cm extended length cuffs for additional coverage of the forearm.</p> <p>Colour: Blue</p> <p>Nominal fingertip thickness: 11mil (.28mm). Must be puncture Resistant.</p> <p>Must have a standard box configuration to allow it to fit into a standard glove box holder.</p> <p>50/Box</p> | | |
| ICN NO: | ITEM: | DESCRIPTION: |
| | NitGlove1 | Non-Sterile Nitrile Decontamination Gloves Size: Small |
| | NitGlove2 | Non-Sterile Nitrile Decontamination Gloves Size: Medium |
| | NitGlove3 | Non-Sterile Nitrile Decontamination Gloves Size: Large |
| | NitGlove4 | Non-Sterile Nitrile Decontamination Gloves Size: X-Large |

STERILE EYE WASH SOLUTION

Purpose: For the irrigation and cleansing of the eyes that have come into contact with chemicals

Consists of a squeeze/pour bottle containing isotonic sterile water.

Must be sealed with a visible expiry date

| ICN NO: | ITEM: | DESCRIPTION: |
|----------------|--------------|---|
| | EyeWash1 | Sterile Eye Wash Solution 200- 250ml |
| | EyeWash2 | Sterile Eye Wash Solution 500ml |

SECTION M

SBD 3.2

PRICING SCHEDULE – NON-FIRM PRICES
(PURCHASES)

NOTE: PRICE ADJUSTMENTS WILL BE ALLOWED AT THE PERIODS AND TIMES SPECIFIED IN THE BIDDING DOCUMENTS.

| | |
|--|------------------------------------|
| Name of bidder..... | Bid number: ZNB 6745/2024-H |
| Closing Time: 11:00 | Closing date: |
| BID DESCRIPTION | |
| DESCRIPTION: HE SUPPLY AND DELIVERY OF HOSPITAL SUNDRIES FOR VARIOUS INSTITUTIONS: 3 YEARS CONTRACT | |

OFFER TO BE VALID FOR 180 DAYS FROM THE CLOSING DATE OF BID.

1. GENERAL SUNDRIES

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|--------------------------------------|-----|------|-----------|--|-----------------------------------|-----------------------------------|-----------------------------------|---|
| PERSONAL IDENTIFICATION BANDS | 1 | | 30 028 07 | Personal Identification Band - Adult Yellow Locking device from 12cm - 25cm | | | | |
| | 2 | | 30 028 08 | Personal Identification Band - Adult Pink Locking device from 12cm - 25cm | | | | |
| | 3 | | 30 028 09 | Personal Identification Band - Adult White Locking device from 12cm - 25cm | | | | |
| | 4 | | 30 028 10 | Personal Identification Band - Child White Locking device from 9cm - 21cm | | | | |
| | 5 | | 30 028 12 | Personal Identification Band - Child Pink Locking device from 9cm - 21cm | | | | |
| | 6 | | 30 028 13 | Personal Identification Band - Child Yellow Locking device from 9cm - 21cm | | | | |
| | 7 | | 30 028 11 | Personal Identification Band - Neonate and Premature Babies White Locking device from 6cm - 12cm | | | | |

| | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|----------------------|-----|------|-----------|--|-----------------------------------|-----------------------------------|-----------------------------------|---|
| DISPOSABLE UNDERWEAR | 8 | | 39 054 13 | Disposable underwear - Small Briefs Waist circumference: 80-90cm Leg circumference: 28-60cm Waist-crotch: ≥ 30cm | | | | |
| | 9 | | 39 054 14 | Disposable underwear - Medium Briefs Waist circumference: 90-100cm Leg circumference: 32-64cm Waist-crotch: ≥ 32cm | | | | |
| | 10 | | 39 054 15 | Disposable underwear - Large Briefs Waist circumference: 100-120cm Leg circumference: 36-72cm Waist-crotch: ≥ 35cm | | | | |
| | 11 | | 39 054 16 | Disposable underwear - X-Large Briefs Waist circumference: 120-140cm Leg circumference: 36-76cm Waist-crotch: ≥ 40cm | | | | |
| | 12 | | 39 054 17 | Disposable underwear - XX-Large Briefs Waist circumference: 130cm-150cm Leg circumference: 40cm-80cm Waist-crotch: 44cm | | | | |
| | 13 | | 39 054 18 | Disposable underwear - XXX-Large Briefs Waist circumference: 140cm-160cm Leg circumference: 44cm-88cm Waist-crotch: 46cm | | | | |
| FLUSHABLE COVERS | 14 | | 39 050 11 | Bedpan Cover - Flushable Size: 30.5cm x 7.6cm x 40.6cm | | | | |
| | 15 | | 39 050 12 | Urinal Covers - Flushable Size: 15.22cm x 7.6cm x 30.5cm | | | | |
| PROTECTIVE COVERS | 16 | | USPcover | Ultrasound probe cover Size: ± 10cm x ± 60cm | | | | |
| | 17 | | 30 461 10 | Temperature probe cover Length: 15-20cm | | | | |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|----------------|-----|------|-----------|---|-----------------------------------|-----------------------------------|-----------------------------------|---|
| WIPES | 18 | | 30 568 79 | Alcohol Swabs - 24mm x 30mm | | | | |
| | 19 | | AdhWipe | Wipe to remove adhesiveness from the skin. Size: 40cm x 40cm | | | | |
| GELS & CREAMS | 20 | | 30 344 30 | Ultrasound Gel with dispenser 200ml | | | | |
| | 21 | | 30 344 32 | Ultrasound Gel Jerry Can 5 Litre | | | | |
| | 22 | | LubGel | Lubricating (KY) Jelly 2.5g Sachets | | | | |
| | 23 | | 53 026 10 | Anaesthetic Gel 2% 20ml Tube | | | | |
| ECG ELECTRODES | 24 | | 30 248 03 | Solid Gel ECG Electrode - Adult | | | | |
| | 25 | | 30 248 04 | Solid Gel ECG Electrode - Paediatric < 50mm wide | | | | |
| | 26 | | 30 248 09 | Sponge ECG Electrode - Paediatric < 50mm wide | | | | |
| | 27 | | 30 248 18 | Sponge ECG Electrode - Adult | | | | |
| | 28 | | 30 248 27 | Tape ECG Electrode - Paediatric < 50mm wide | | | | |
| | 29 | | 30 248 37 | Tape ECG Electrode - Neonate < 30mm wide | | | | |
| | 30 | | 30 248 38 | Tape ECG Electrode - Adult | | | | |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|-----------------|-----|------|-----------|---|-----------------------------------|-----------------------------------|-----------------------------------|---|
| PLASTIC APRONS | 31 | | 39 010 15 | Full Body Apron Medium - General use Length from neck: ± 110cm Colour: Blue | | | | |
| | 32 | | 39 010 16 | Full Body Apron Large - General use Length from neck: ± 130cm Colour: Blue | | | | |
| | 33 | | 39 010 18 | Full Body Apron Medium - Sluice room Length from neck: ± 110cm Colour: Red | | | | |
| | 34 | | 39 010 19 | Full Body Apron Large - Sluice room Length from neck: ± 130cm Colour: Red | | | | |
| | 35 | | 39 010 21 | Full Body Apron Medium - Infectious areas Length from neck: ± 110cm Colour: Yellow | | | | |
| | 36 | | 39 010 22 | Full Body Apron Large - Infectious areas Length from neck: ± 130cm Colour: Yellow | | | | |
| | 37 | | 39 010 24 | Full Body Apron Medium - Clean areas (kitchen) Length from neck: ± 110cm Colour: White | | | | |
| | 38 | | 39 010 25 | Full Body Apron Large - Clean areas (kitchen) Length from neck: ± 130cm Colour: White | | | | |
| | 39 | | 39 010 26 | Surgical Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Light Green | | | | |
| | 40 | | 54 020 10 | Laundry Heavy Duty Full Body Apron - Reusable Protective length: ± 120cm Colour: Red | | | | |
| NON-DISPOSABLES | 41 | | 30 542 06 | Stethoscope – Infant Diaphragm: ± 2.7cm | | | | |
| | 42 | | 30 542 07 | Stethoscope – Paediatric Diaphragm: ± 3.3cm | | | | |
| | 43 | | 30 542 08 | Stethoscope – Adult Diaphragm: ± 5cm | | | | |
| | 44 | | 30 244 00 | Thermometer – Digital temperature ranges of at least: 34°C to 42°C. | | | | |
| | 45 | | ND-PB | Pressure Bag – Non-Disposable | | | | |
| | 46 | | ET-M | Manometer - Endotracheal and Supraglottic Cuff Pressure Monitor | | | | |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|----------------|-----|------|-----------|---|-----------------------------------|-----------------------------------|-----------------------------------|---|
| PASTEURISATION | 47 | | 30 047 00 | Pasteurisation Bottle - Glass. Volume: 125ml | | | | |
| | 48 | | 30 047 01 | Lid for Glass Pasteurisation Bottle 48mm single thread | | | | |
| | 49 | | 30 047 02 | Pasteurisation and Human Milk Storage Bottle - BPA-free Plastic Volume: 130ml ; Neck Diameter: 50mm | | | | |
| | 50 | | 30 047 03 | Lids with foil seal for Human Milk Storage Bottles - Reusable | | | | |
| | 51 | | 30 052 09 | Bottle Brush: Length including brush: 40cm . Brush: 10cm , Handle: 30cm Flexible bristles: 5cm | | | | |
| | 52 | | PL1 | Small-scale Pasteuriser Labels Size: 28mm x 89mm | | | | |
| | 53 | | PL2 | Plasticised Pasteuriser Labels Size: 100mm x 25mm | | | | |
| NEONATAL | 54 | | 30 445 23 | Litmus Paper Dark Blue | | | | |
| | 55 | | 30 445 24 | Litmus Paper Red | | | | |
| | 56 | | 30 445 25 | pH Indicator Test strips | | | | |
| | 57 | | 30 460 51 | Phototherapy Eye-Shield - Preemie , head circumference ± 26-32cm | | | | |
| | 58 | | 30 460 52 | Phototherapy Eye-Shield - Regular , head circumference ± 33-38 cm | | | | |
| | 59 | | 30 072 01 | Umbilical Cord Clamp. ± 5cm long | | | | |
| | 60 | | 30 326 10 | Amniotomy hook - 25cm | | | | |
| CONNECTORS | 61 | | 30 089 23 | Y-Shaped Tubing Connector - Large OD: ±10mm All sides must be of equal length – 30mm length from point to junction | | | | |
| | 62 | | 30 089 36 | Straight Tubing Connector - Open-ended. Double tapered Small OD: 7-11mm; ID: 5mm; Total length: 6cm | | | | |
| | 63 | | 30 089 37 | Straight Tubing Connector - Open-ended. Double tapered Large OD: 10 - 14mm ID: 8mm Total length: 6cm | | | | |
| | 64 | | 30 534 80 | Spigot - Small Diameter: 9mm; Length: ±40mm | | | | |
| | 65 | | 30 534 82 | Spigot - Large Diameter: 10mm; Length: ± 52mm | | | | |

| CATEGORY: | NO: | ICN: | ITEM: | DESCRIPTION: | Unit PRICE (INCL VAT) for year 1 | Unit PRICE (INCL VAT) for year 2 | Unit PRICE (INCL VAT) for year 3 | Sub-Total Price =SUM OF YR 1,2 &3 (incl. VAT) |
|--|-----|------|-----------|--|-----------------------------------|-----------------------------------|-----------------------------------|---|
| DRESSING PACKS | 66 | | 30 432 15 | Large Tray Dressing Pack for Minor Surgical Procedures | | | | |
| | 67 | | 30 432 16 | Multi-Pack Tray Dressing Pack for Minor Surgical Procedures | | | | |
| | 68 | | 30 432 17 | Catheter Tray Dressing Pack | | | | |
| | 69 | | 30 432 18 | Maternity Swabbing Tray Dressing Pack | | | | |
| COTTON WOOL AND CLEANSING SWABS | 70 | | 30 645 05 | Cotton wool roll – 1st Grade 50g roll | | | | |
| | 71 | | 30 645 06 | Cotton wool roll – 2nd Grade 50g roll | | | | |
| | 72 | | 30 568 03 | Cotton Wool Absorbent Balls – Sterile 5 balls Sterile in an easy to open peel pouch | | | | |
| | 73 | | 30 568 04 | Cotton Wool Absorbent Balls – Unsterile Box of 1000 balls | | | | |
| | 74 | | 30 568 05 | Cleansing ball swabs - X-Ray detectable, Unsterile, Diameter: 65mm | | | | |
| | 75 | | 30 568 06 | Cleansing ball swabs - X-Ray detectable, Sterile; Diameter: 65mm | | | | |

- Required by: **VARIOUS INSTITUTIONS**
- Country of origin:
- Delivery Period:

Note: All delivery costs must be included in the bid price, for delivery at the prescribed destination.

.....
 (Signature of Bidder) Date (Signature of Witness) Date