



Daily Equipment Checklist-General. Reg./Tert.

Month and Year:

Unit:

Check each piece of equipment daily.

- Cleaned daily
- Fully Functional-no error readings
- All consumables eg probes/cables present and functional
- No cracks/breaks
- No exposed wires
- Phototherapy lights <1000hrs & all working

Equipment item:	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Breast pump																															
Cold light																															
Diagnostic set																															
Double O ₂ flowmeter																															
Flash bilirubinometer																															
Fridge-Immunisation																															
Fridge-Medication																															
Fridge-Milk																															
Glucometer																															
Haemoglobinometer																															
Neo-puff																															
Phototherapy blanket																															
Phototherapy lights																															
Portable O ₂ cylinder with pressure gauge																															
Portable suction																															
Scale infant electronic																															
Stethoscope																															
Transport incubator																															
Wall thermometer																															
X-ray viewing box – long																															
X-ray viewing box – short																															
12 lead ECG machine																															
Blood gas machine																															
Blood warmer																															
Head/body cooling system																															
Lactate meter																															
Oto acoustic emission																															
Ret-cam																															
US machine																															
Cardiac output monitor																															
Nitric oxide																															

